



Innovation. Protection. Comfort.

Secure-Gard[®] and Insta-Gard[®] Facial Protection

BY CONVERTORS/SURGICAL SPECIALTIES

 **Allegiance**

a Cardinal Health company

WHEN YOU CHOOSE ALLEGIANCE HEALTHCARE FOR YOUR MEDICAL SUPPLY NEEDS, YOU KNOW YOU ARE GETTING TOP-QUALITY PRODUCTS THAT OFFER OUTSTANDING PERFORMANCE AND VALUE. OUR ABILITY TO OFFER INNOVATIVE PRODUCTS COMES FROM YEARS OF WORKING TO MEET THE NEEDS OF HEALTHCARE PROFESSIONALS AROUND THE WORLD.

THAT'S WHY WE ARE PROUD TO INTRODUCE THE SECURE-GARD® AND INSTA-GARD® FACE MASKS, FACE SHIELDS, AND VISOR SYSTEMS. FROM PRIMARY CARE TO SURGICAL SETTINGS, YOU CAN WORK CONFIDENTLY, KNOWING THAT YOU'VE CHOSEN A PRODUCT THAT WILL HELP PROTECT YOU AND YOUR STAFF FROM A WIDE RANGE OF PATHOGENS.

EACH SECURE-GARD AND INSTA-GARD FACE MASK IS MANUFACTURED IN A CERTIFIED ISO 9001 FDA REGISTERED FACILITY, AND UNDERGOES RIGOROUS QUALITY INSPECTIONS AT EVERY STEP OF THE MANUFACTURING PROCESS TO ENSURE THAT THE PRODUCT YOU RECEIVE IS THE PRODUCT YOU EXPECT — THE PRODUCT YOU DESERVE.

DISCOVER FOR YOURSELF THE QUALITY AND WORKMANSHIP THAT SETS SECURE-GARD AND INSTA-GARD FACIAL PROTECTION APART FROM OTHER PRODUCTS IN THE MARKET. FOR A FREE SAMPLE, CONTACT YOUR ALLEGIANCE HEALTHCARE SALES REPRESENTATIVE TODAY, OR CALL US DIRECTLY AT [1.800.964.5227](tel:18009645227).

Innovation and Technology.

- Proprietary eye shield attachment positions the shield away from the user's face for a more comfortable fit
- Pre-scored (folded) eye shield provides better wraparound face coverage

Protection.

- Downward-folding pleats help minimize splashed fluids from penetrating the mask
- Engineered materials to deliver a high level of filtration efficiency and fluid resistance protection

Comfort.

- Three-layer construction allows for better breathability than most four-layer masks, without sacrificing filtration efficiency
- Contoured nose cutout on all eye shields enhances breathing comfort



FILTRATION



FOG-FREE



FLUID RESISTANT



EYE SHIELD

Choose the mask that best fits your applications

Earloop

WITH FILTRATION

<input type="checkbox"/> WITH FOG-FREE	<input checked="" type="checkbox"/> WITH FLUID RESISTANT	<input checked="" type="checkbox"/> WITH EYE SHIELD Cat. # AT74631 (pg. 14)
	<input type="checkbox"/> WITHOUT FLUID RESISTANT	<input type="checkbox"/> WITHOUT EYE SHIELD Cat. # AT74531 (pg. 14)
<input type="checkbox"/> WITHOUT FOG-FREE	<input checked="" type="checkbox"/> WITH FLUID RESISTANT	<input checked="" type="checkbox"/> WITH EYE SHIELD Cat. # AT71311 (pg. 14)
	<input type="checkbox"/> WITHOUT FLUID RESISTANT	<input type="checkbox"/> WITHOUT EYE SHIELD Cat. # AT7511-FF (pg. 14)
<input type="checkbox"/> WITHOUT FOG-FREE	<input checked="" type="checkbox"/> WITH FLUID RESISTANT	<input checked="" type="checkbox"/> WITH EYE SHIELD Cat. # AT7515-WE (pg. 14)
	<input type="checkbox"/> WITHOUT FLUID RESISTANT	<input type="checkbox"/> WITHOUT EYE SHIELD Cat. # AT7515 (pg. 14)
<input type="checkbox"/> WITHOUT FOG-FREE	<input checked="" type="checkbox"/> WITH FLUID RESISTANT	<input checked="" type="checkbox"/> WITH EYE SHIELD Cat. # AT7511-WE (pg. 14)
	<input type="checkbox"/> WITHOUT FLUID RESISTANT	<input type="checkbox"/> WITHOUT EYE SHIELD Cat. # AT7511-W (pg. 11), AT70021 (pg. 13), AT71021 (pg. 13)

Surgical Tie

WITH FILTRATION

<input checked="" type="checkbox"/> WITH FOG-FREE	<input checked="" type="checkbox"/> WITH FLUID RESISTANT	<input checked="" type="checkbox"/> WITH EYE SHIELD Cat. # AT74935 (pg. 8), AT74635 (pg. 10)
	<input type="checkbox"/> WITHOUT FLUID RESISTANT	<input type="checkbox"/> WITHOUT EYE SHIELD Cat. # AT744535 (pg. 8)
<input type="checkbox"/> WITHOUT FOG-FREE	<input checked="" type="checkbox"/> WITH FLUID RESISTANT	<input checked="" type="checkbox"/> WITH EYE SHIELD Cat. # AT71335 (pg. 10)
	<input type="checkbox"/> WITHOUT FLUID RESISTANT	<input type="checkbox"/> WITHOUT EYE SHIELD Cat. # AT72835 (pg. 8), AT71235 (pg. 10)
<input type="checkbox"/> WITHOUT FOG-FREE	<input checked="" type="checkbox"/> WITH FLUID RESISTANT	<input checked="" type="checkbox"/> WITH EYE SHIELD Cat. # AT7506-WE (pg. 10)
	<input type="checkbox"/> WITHOUT FLUID RESISTANT	<input type="checkbox"/> WITHOUT EYE SHIELD Cat. # AT7506 (pg. 10)
<input type="checkbox"/> WITHOUT FOG-FREE	<input checked="" type="checkbox"/> WITH FLUID RESISTANT	<input checked="" type="checkbox"/> WITH EYE SHIELD Cat. # AT71135 (pg. 10)
	<input type="checkbox"/> WITHOUT FLUID RESISTANT	<input type="checkbox"/> WITHOUT EYE SHIELD Cat. # AT71035 (pg. 7), AT7505-P (pg. 9)

Glossary

ADHESIVE TAPE FOG-FREE STRIP

Double-sided adhesive tape designed to temporarily adhere to the dermal layer, in an effort to minimize warm exhaled air from occluding eye wear or eye shields.

ANTI-GLARE STRIP

The black (or dark-colored) strip placed on the top of a face mask across the nose section to help absorb and redirect bright light away from the wearer's eyes.

BFE % (BACTERIAL FILTRATION EFFICIENCY)

Measures the percent efficiency at which the face mask stops bacteria from passing through the mask.

BLUE FOAM FOG-FREE STRIP

Open cell polyester that is designed to provide a foam barrier between the top inner layer of the face mask and the bridge of the nose and cheekbones, in an effort to minimize warm exhaled air from occluding eye wear or eye shields.

DELTA P (ΔP)

The pressure drop across a face mask, expressed in mm water/cm² — the higher the Delta P, the more difficult the mask is to breathe through.

FLAT WHITE FOG-FREE STRIP

Meltblown polypropylene designed to lay flat and bridge the gap between the top of the mask and the top inner layer of the face mask and the bridge of the nose and cheekbones, in an effort to minimize warm exhaled air from occluding eye wear or eye shields.

FLUID REPELLENT

The ability of a face mask's material construction to minimize fluids from traveling through the material and potentially coming into contact with the user of the face mask. Material constructions that are fluid repellent typically do not minimize the transfer of fluids through the material as well as Fluid Resistant constructions.

FLUID RESISTANT

The ability of a face mask's material construction to minimize fluids from traveling through the material and potentially coming into contact with the user of the face mask.

ISOLATION MASK

A general purpose mask that is usually donned when transporting patients or soiled linens that have not been designated high-risk. Typically, isolation masks offer a minimum level of protection to the wearer.

PFE % (PARTICULATE FILTRATION EFFICIENCY)

Measures the percent efficiency at which the face mask stops particulate matter from passing through the mask.

PROCEDURE MASK

Another term for an earloop face mask. These masks are designed to be used primarily in the healthcare setting for general procedures and checkups. They are designed to be changed in between each patient contact.

SONICALLY SEALED

A bond on a face mask that is produced by utilizing ultrasonic sound waves to form a seal.

SURGICAL MASK

Another term for a healthcare face mask with ties as opposed to earloops. These masks are designed to be used primarily in the healthcare setting during standard surgical procedures. Since surgical masks are worn with ties as opposed to earloops, they generally offer a more secure and personalized fit for the wearer; minimizing the potential of bacterial or particulate penetration.

TORTUOUS PATHWAY

This pathway can be created by placing several layers of nonwoven material together. This construction makes it difficult for bacteria and particles to travel through the material.

WHITE CUSHIONED FOG-FREE STRIP

Closed cell polyethylene designed to offer a smooth cushioned barrier between the top inner layer of the face mask and the bridge of the nose and cheekbones, in an effort to minimize warm exhaled air from occluding eye wear or eye shields.



FROM THE OPERATING ROOM TO THE RESEARCH LAB, ALLEGIANCE HEALTHCARE'S SECURE-GARD® AND INSTA-GARD® FACE MASKS CAN HELP PROTECT YOU FROM A WIDE RANGE OF PATHOGENS. AVAILABLE IN AN ARRAY OF STYLES AND COLORS, WE HAVE A MASK TO MEET YOUR SPECIFIC NEEDS.



Secure-Gard[®] Surgical Masks with Surgical Ties

SECURE-GARD[®] STANDARD SURGICAL MASKS

Secure-Gard surgical masks feature ties and bindings that are sonically sealed as opposed to being sewn or glued, which helps maximize both comfort and durability. These masks are an excellent choice for non-invasive surgical procedures where there is little or no possibility of fluids being splashed or sprayed. Extra long ties provide a more secure fit, which helps to minimize "blow by."

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT71035 50 per box, 300 per case	Blue vertical tie surgical mask with soft, breathable inner facing	≥97%	≥98%	≥98%	≤2.0	N/A



FILTRATION



FOG-FREE



FLUID RESISTANT



EYE SHIELD

SECURE-GARD® SURGICAL MASKS WITH FOG-FREE STRIP

Minimize fogging and improve your level of clarity with one of the Secure-Gard surgical masks with fog-free strip. Available with adhesive strip or white cushioned foam strip, these masks provide an almost impenetrable, contoured seal between the mask and your face to minimize warm air (breath) coming in contact with glasses, microscopes, or protective eyewear. Adhesive strips adhere to the face and contour around the nose and cheekbones for a customized fit for everyone in the facility.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT741435 50 per box, 300 per case	Blue vertical tie surgical mask with white cushioned fog-free strip	≥97%	≥98%	≥98%	≤2.0	N/A
AT72835 50 per box, 300 per case	Green vertical tie surgical mask with adhesive tape fog-free strip	≥97%	≥98%	≥98%	≤2.0	N/A

SECURE-GARD® FLUID RESISTANT SURGICAL MASKS WITH FOG-FREE STRIP

Designed to offer a high level of facial protection, the Secure-Gard surgical mask with fluid resistance exceeds ASTM standard F1862 for splash resistance.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT744535 50 per box, 300 per case	Apricot fluid resistant , vertical tie surgical mask with white cushioned fog-free strip	≥99%	≥99%	≥99%	≤4.8	160

SECURE-GARD® FLUID RESISTANT SURGICAL MASKS WITH FOG-FREE STRIP AND WRAPAROUND EYE SHIELD

Protect your eyes, nose, and mouth in one easy step with a specially designed Secure-Gard facemask with wraparound eye shield. Enhanced with our proprietary inner-media layer that minimizes fluid penetration while filtering out bacteria. All shields are treated with an anti-fogging agent and are pre-scored for excellent coverage.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT744635 25 per box, 100 per case	Apricot fluid resistant, vertical tie surgical mask with wraparound eye shield, white cushioned fog-free strip , and broad anti-glare strip	≥99%	≥99%	≥99%	≤4.8	160
AT74935 25 per box, 100 per case	Apricot fluid resistant, vertical tie surgical mask with adhesive tape fog-free strip , wraparound eye shield, and broad anti-glare strip	≥99%	≥99%	≥99%	≤4.8	160



Insta-Gard[®] Surgical Masks with Surgical Ties

INSTA-GARD[®] STANDARD SURGICAL MASKS

Designed to offer a soft, comfortable fit, Insta-Gard surgical tie masks offer the protection you need at a price that can fit into any budget. Designed with downward folding pleats and an enclosed nosepiece, Insta-Gard surgical masks are an excellent choice for minimally invasive surgical procedures.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT7505-P 50 per box, 300 per case	Blue vertical tie surgical mask	≥96%	≥97%	≥97%	≤2.0	N/A
AT7505-PG 50 per box, 300 per case	Green vertical tie surgical mask	≥96%	≥97%	≥97%	≤2.0	N/A



FILTRATION



FOG-FREE



FLUID RESISTANT



EYE SHIELD

INSTA-GARD® SURGICAL MASKS WITH FOG-FREE STRIP

If you like the standard Insta-Gard surgical tie face mask, you will appreciate the standard mask with a foam fog-free strip. Available with different inner facings and two different types of fog-free strips to meet a variety of needs, these masks provide protection while helping to minimize fogging.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT7505-PFF 50 per box, 300 per case	Blue vertical tie surgical mask with flat white fog-free strip and soft, breathable inner facing	≥96%	≥97%	≥97%	≤2.0	N/A
AT71039 50 per box, 600 per case	Blue horizontal tie surgical mask with flat white fog-free strip	≥96%	≥97%	≥97%	≤2.0	N/A
AT71235 50 per box, 300 per case	Blue vertical tie surgical mask with blue foam fog-free strip and soft, breathable inner facing	≥96%	≥97%	≥97%	≤2.0	N/A
AT72235 50 per box, 300 per case	Green vertical tie surgical mask with blue foam fog-free strip	≥96%	≥97%	≥97%	≤2.0	N/A

INSTA-GARD® FLUID RESISTANT SURGICAL MASKS WITH FOG-FREE STRIP

Insta-Gard surgical masks offer 120mm Hg ASTM-level fluid resistance protection while providing excellent filtration efficiencies. Available with or without an eye shield, these masks help protect your staff and your budget at the same time.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT74535 50 per box, 300 per case	Apricot fluid resistant , vertical tie surgical mask with blue foam fog-free strip	≥98%	≥98%	≥98%	≤3.6	120

INSTA-GARD® FLUID RESISTANT SURGICAL MASKS WITH FOG-FREE STRIP AND WRAPAROUND EYE SHIELD

These Insta-Gard face masks feature an eye shield pre-scored (folded) on the sides to offer excellent protection for the side of the face, and a contoured nose cut-out to help exhaled breath escape unobstructed through the mask. For great protection at a great price, try an Insta-Gard face mask with eye shield.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT74635 25 per box, 100 per case	Apricot fluid resistant , vertical tie surgical mask with blue foam fog-free strip , broad anti-glare strip, and wraparound eye shield	≥98%	≥98%	≥98%	≤3.6	120

ADDITIONAL INSTA-GARD® SURGICAL MASKS

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT71135 25 per box, 100 per case	Blue vertical tie surgical mask with wraparound eye shield	≥96%	≥97%	≥97%	≤2.0	N/A
AT71335 25 per box, 100 per case	Blue vertical tie surgical mask with blue foam fog-free strip , wraparound eye shield , and anti-glare strip	≥96%	≥97%	≥97%	≤2.0	N/A
AT7506 50 per box, 300 per case	Teal fluid resistant , vertical tie surgical mask	≥98%	≥98%	≥98%	≤3.6	120
AT7506-WE 25 per box, 100 per case	Teal fluid resistant , vertical tie surgical mask with wraparound eye shield	≥98%	≥98%	≥98%	≤3.6	120



Secure-Gard[®] Procedure Masks with Earloops

SECURE-GARD[®] STANDARD PROCEDURE MASKS

Secure-Gard procedure masks offer excellent filtration efficiency and feature earloops for quick and easy donning and doffing. Secure-Gard procedure masks are a great choice for patient care procedures and the ER when there is no expectation of fluids being splashed or sprayed. Available in several colors.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT7511 50 per box, 300 per case	Blue procedure mask with earloops	≥97%	≥98%	≥98%	≤2.0	N/A
AT7511-PK 50 per box, 300 per case	Pink procedure mask with earloops	≥97%	≥98%	≥98%	≤2.0	N/A
AT7511-G 50 per box, 300 per case	Green procedure mask with earloops	≥97%	≥98%	≥98%	≤2.0	N/A
AT7511-W 50 per box, 300 per case	White procedure mask with earloops	≥97%	≥98%	≥98%	≤2.0	N/A



FILTRATION



FOG-FREE



FLUID RESISTANT



EYE SHIELD

SECURE-GARD® PROCEDURE MASKS WITH FOG-FREE STRIP

Experience the comfort and clarity of Secure-Gard procedure masks with fog-free strip. Earloops offer quick and easy donning and doffing.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT741411 50 per box, 300 per case	Blue procedure mask with white cushioned fog-free strip	≥97%	≥98%	≥98%	≤2.0	N/A

SECURE-GARD® FLUID RESISTANT PROCEDURE MASKS WITH FOG-FREE STRIP

Don't settle for uncomfortable polyester foam strips when you can get the fluid resistance you need in a mask with a soft, comfortable, cushioned fog-free strip. Experience the comfort and dependability of Secure-Gard today.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT744531 25 per box, 100 per case	Apricot fluid resistant procedure mask with white cushioned fog-free strip	≥99%	≥99%	≥99%	≤4.8	I60

SECURE-GARD® FLUID RESISTANT PROCEDURE MASKS WITH FOG-FREE STRIP AND WRAPAROUND EYE SHIELD

Secure-Gard procedure masks with wraparound eye shield offer soft, comfortable earloops and three-layer construction for hours of comfortable wear. All masks feature a contoured nose cutout for ease of breathing and a broad anti-glare strip.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT744631 25 per box, 100 per case	Apricot fluid resistant procedure mask with white cushioned fog-free strip , wraparound eye shield, and broad anti-glare strip	≥99%	≥99%	≥99%	≤4.8	I60

FOR ADDITIONAL INFORMATION, PLEASE CONTACT YOUR ALLEGIANCE HEALTHCARE SALES REPRESENTATIVE OR CUSTOMER SERVICE AT [1.800.964.5227](tel:1.800.964.5227).

*microns



Insta-Gard[®] Procedure Masks with Earloops

INSTA-GARD[®] STANDARD PROCEDURE MASKS

Designed to offer more protection than a standard isolation mask, Insta-Gard procedure masks are available in blue or yellow with various inner facings to fit your personal preference.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT71021 50 per box, 500 per case	Blue procedure mask with earloops	≥96%	≥97%	≥97%	≤2.0	N/A
AT70021 50 per box, 500 per case	Yellow procedure mask with earloops	≥96%	≥97%	≥97%	≤2.0	N/A
AT7311 50 per box, 300 per case	Blue procedure mask with earloops	≥96%	≥97%	≥97%	≤2.0	N/A
AT7311-Y 50 per box, 300 per case	Yellow procedure mask with earloops	≥96%	≥97%	≥97%	≤2.0	N/A



FILTRATION



FOG-FREE



FLUID RESISTANT



EYE SHIELD

INSTA-GARD® PROCEDURE MASKS WITH FOG-FREE STRIP

Available with standard blue foam fog-free strip, Insta-Gard procedure masks offer great protection for you and your staff, at an affordable price. Available with or without an eye shield, this mask is a perfect choice for quick procedures and patient exams.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT7511-FF 50 per box, 300 per case	Blue procedure mask with blue foam fog-free strip	≥96%	≥97%	≥97%	≤2.0	N/A

INSTA-GARD® FLUID RESISTANT PROCEDURE MASKS WITH FOG-FREE STRIP

All fluid resistant Insta-Gard masks meet the ASTM standard of protection at the 120mm Hg level, and are comfortable during extended hours of wear. Insta-Gard procedure masks are available with or without an eye shield and in various colors to fit your needs.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT74531 50 per box, 300 per case	Apricot fluid resistant procedure mask with blue foam fog-free strip	≥98%	≥98%	≥98%	≤3.6	120

INSTA-GARD® FLUID RESISTANT PROCEDURE MASKS WITH FOG-FREE STRIP AND WRAPAROUND EYE SHIELD

Insta-Gard procedure masks with wraparound eye shield offer sturdy construction and excellent breathability for hours of comfortable wear. All masks feature a contoured nose cutout for ease of breathing and a broad anti-glare strip.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT74631 25 per box, 100 per case	Apricot fluid resistant procedure mask with blue foam fog-free strip and wraparound eye shield with broad anti-glare strip	≥98%	≥98%	≥98%	≤3.6	120

ADDITIONAL INSTA-GARD® PROCEDURE MASKS

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT71311 25 per box, 100 per case	Blue procedure mask with blue foam fog-free strip and wraparound eye shield	≥96%	≥97%	≥97%	≤2.0	120
AT7511-WE 25 per box, 100 per case	Blue procedure mask with wraparound eye shield	≥96%	≥97%	≥97%	≤2.0	N/A
AT7515 50 per box, 300 per case	Teal fluid resistant procedure mask	≥98%	≥98%	≥98%	≤3.6	120
AT7515-WE 25 per box, 100 per case	Teal fluid resistant procedure mask with wraparound eye shield	≥98%	≥98%	≥98%	≤3.6	120



Secure-Gard[®] Specialty Masks

SECURE-GARD[®] CONE MASKS

Shaped to fit securely around the nose and mouth, Secure-Gard cone masks offer superb protection from a wide range of bacteria while maintaining comfort and enhanced breathability.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT7509 50 per box, 300 per case	Fluid resistant cone mask, molded cone shape	≥99.5% ¹	N/A	N/A	≤0.3	N/A

SECURE-GARD[®] LASER MASKS

The Secure-Gard laser face mask is designed to fit snugly around the nose and mouth while offering excellent protection against inhalation of small particles introduced into the air during laser surgery. Don't let the soft feel of this mask fool you. Secure-Gard laser masks work hard to protect you. Available in Cool Grey color.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT75435 50 per box, 300 per case	Grey fluid resistant vertical tie laser surgery mask	≥99%	≥99%	≥99%	≤4.8	160



FILTRATION



FOG-FREE



FLUID RESISTANT



EYE SHIELD

* microns ¹ Test based on Green & Vesley Test Method results.

SECURE-GARD® N95 PARTICULATE RESPIRATOR

Manufactured in compliance with NIOSH Standards for N95 respirators, the Secure-Gard N95 respirator is soft and comfortable while providing a secure fit around your nose and mouth. Can be used in settings where TB is believed to be present or during surgery. This mask offers excellent protection to the wearer.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.3* ²	Delta P	ASTM Level
AT79595 20 per box, 120 per case	Type N95 NIOSH Certified Particulate Respirator	—	—	≥95%	≤12	—

SECURE-GARD® BRAND SECURE BILL™ SURGICAL MASKS

Minimize that close, closed-in feeling with one of our Secure Bill surgical tie masks. Designed to stand away from the nose and mouth, these masks provide the wearer a more open breathing sensation without compromising protection. Standard and fluid resistant Secure Bill masks are both available with or without eye shields.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT51035 50 per box, 300 per case	Blue Secure Bill surgical mask with soft, breathable inner facing	≥97%	≥98%	≥98%	≤2.0	N/A
AT51135 25 per box, 100 per case	Blue Secure Bill surgical mask with soft, breathable inner facing and wraparound eye shield	≥97%	≥98%	≥98%	≤2.0	N/A
AT54535 50 per box, 300 per case	Apricot fluid resistant Secure Bill surgical mask with white cushioned fog-free strip	≥99%	≥99%	≥99%	≤4.8	160
AT54635 25 per box, 100 per case	Apricot fluid resistant Secure Bill surgical mask with white cushioned fog-free strip and wraparound eye shield	≥99%	≥99%	≥99%	≤4.8	160

SECURE-GARD® SENSITIVE SKIN SURGICAL MASK

Manufactured utilizing a unique blend of materials, the Secure-Gard Sensitive Skin mask is designed to be one of the softest, most comfortable masks available in the healthcare marketplace today. The Secure-Gard Sensitive Skin surgical mask combines excellent protection and comfort to provide you the very best a face mask has to offer.

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.3*	Delta P	ASTM Level
AT73035 50 per box, 300 per case	White sensitive skin surgical tie mask	≥97%	≥98%	≥98%	≤2.0	N/A

*microns ²Test results based on NACL for Particulates.

Highly Breathable/High Filtration Face Masks

Manufactured utilizing a state-of-the-art material that provides excellent breathability without sacrificing protection, our highly breathable/high filtration masks are the perfect choice for long procedures when a mask is absolutely critical. You will be amazed at how easy it is to breathe through this mask. Call your Allegiance Healthcare sales representative for more information or samples.

• All masks designated "Fluid Resistant" meet or exceed ASTM F1862 for fluid resistance

• Unique material construction offers one of the highest levels of breathability and comfort for face masks on the market today

• **Latex-free construction**

• **Meets or exceeds industry standards** (ASTM F23.40.408) for minimizing bacterial filtration

• Features enclosed nose piece for a comfortable, secure fit

• Manufactured in accordance with the Quality System Regulations

STANDARD SURGICAL MASKS WITH TIES

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.3*	Delta P	ASTM Level
AT81075 50 per box, 300 per case	Blue standard face mask with vertical ties	≥98%	≥97%	≥97%	≤0.5	N/A

STANDARD SURGICAL TIE MASKS WITH FOG-FREE STRIP

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.3*	Delta P	ASTM Level
AT81375 40 per box, 240 per case	Blue standard face mask with blue foam fog-free strip	≥98%	≥97%	≥97%	≤0.5	N/A

FLUID RESISTANT SURGICAL TIE MASK WITH FOG-FREE STRIP

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT81335 40 per box, 240 per case	Blue fluid resistant with blue foam fog-free strip	≥98%	≥98%	≥98%	≤1.5	160
AT81535 40 per box, 240 per case	Teal fluid resistant with adhesive tape fog-free strip	≥98%	≥98%	≥98%	≤1.5	160

FLUID RESISTANT SURGICAL TIE MASK WITH FOG-FREE STRIP AND WRAPAROUND EYE SHIELD

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.3*	Delta P	ASTM Level
AT81435 25 per box, 100 per case	Teal fluid resistant with blue foam fog-free strip and wraparound eye shield	≥98%	≥98%	≥98%	≤1.5	160



FILTRATION



FOG-FREE



FLUID RESISTANT



EYE SHIELD

Highly Breathable/High Filtration Procedure Earloop Masks

STANDARD PROCEDURE MASKS WITH EARLOOPS

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
ATH81073 100 per bag, 300 per case	Blue standard procedure mask	≥98%	≥97%	≥97%	≤0.5	N/A

STANDARD PROCEDURE MASKS WITH EARLOOPS AND FOG-FREE STRIP

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT81371 30 per box, 300 per case	Blue standard procedure mask with blue foam fog-free strip	≥98%	≥97%	≥97%	≤0.5	N/A

FLUID RESISTANT PROCEDURE MASKS WITH EARLOOPS AND FOG-FREE STRIP

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT81031 30 per box, 300 per case	Teal fluid resistant procedure mask with blue foam fog-free strip	≥98%	≥98%	≥98%	≤1.5	160

FLUID RESISTANT PROCEDURE MASKS WITH EARLOOPS, FOG-FREE STRIP, AND WRAPAROUND EYE SHIELD

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT81431 25 per box, 100 per case	Teal fluid resistant procedure mask with blue foam fog-free strip and wraparound eye shield	≥98%	≥98%	≥98%	≤1.5	160

Insta-Gard® Isolation Masks with Earloops

INSTA-GARD® STANDARD ISOLATION MASKS

Item Number	Product Description	Filtration % @ 3.0*	Filtration % @ 1.0*	Filtration % @ 0.1*	Delta P	ASTM Level
AT76061 50 per box, 500 per case	Pink isolation mask	≥96%	≥97%	≥97%	≤2.2	N/A
AT70061 50 per box, 500 per case	Yellow isolation mask	≥96%	≥97%	≥97%	≤2.2	N/A
AT71061 50 per box, 500 per case	Blue isolation mask	≥96%	≥97%	≥97%	≤2.2	N/A

*microns

OSHA (Occupational Safety & Health Administration) requires the employer to evaluate the task and type of potential exposure anticipated (level of risk), and to select the appropriate protective apparel that prevents contamination of an employee's skin or clothing from infectious materials.

Cross-Reference

Product Description	Allegiance	Kimberly-Clark/Technol	Precept	3M	Medline	Johnson & Johnson
Standard Procedure Mask	AT7511	47080		N/A	NON27375	4228
Standard Surgical Tie Mask	AT71035 AT71235 AT7505-P	48201	65-3010 65-3100 65-3110 65-3112	N/A	NON27376	4230
Fluid Resistant Surgical Tie Mask	AT74535 AT7506	48247	65-3312	1818	NON27412	4234
Fluid Resistant Surgical Tie Mask with Eye Shield	AT74635 AT74935	48207	N/A	1818FS	NON27410	4232
Fluid Resistant Procedure Mask	AT74531	47108	65-3347	1820	NON27412EL	4224
Fluid Resistant Procedure Mask with Eye Shield	AT74631	47147	65-3345	1820FS	NON27410EL	4267
Laser Mask	AT75435	47625 47650	65-3310 65-3311	1812*	NON24374	4236
N95 Particulate Respirator	AT79595	62126 62355	65-3395	1860 1860S	NON24501 NON24500	4270 4271
Surgical Tie Mask with Adhesive Tape Fog-free Strip	AT72835	49215 49210	65-3322 65-3321	N/A	NON27379A	N/A
Secure Bill™ Surgical Tie Mask (or similar design)	AT51035	375051 48220 375101 48225 375251 49210	65-3124 65-3126	1838	NON27382	4233
Secure Bill™ Surgical Tie Mask with Eye Shield (or similar design)	AT54635	N/A	65-3346	1838FS	N/A	N/A
Isolation Mask	AT71061	47117 47080	65-3450	N/A	NON27120	N/A
Sensitive Skin Mask	AT73035	48390	65-3130	N/A	NON27400	N/A
Cone Mask	AT7509	47200	—	1880+*	NON27381	4235

* This mask is designated Fluid Resistant



a Cardinal Health company

FOR ADDITIONAL INFORMATION, PLEASE CONTACT YOUR ALLEGIANCE HEALTHCARE SALES REPRESENTATIVE OR CUSTOMER SERVICE AT **1.800.964.5227**.

© Copyright 2002, Allegiance Healthcare Corporation.
All rights reserved.
Lit. No. S5000100 (0202/7.5M)

Allegiance Healthcare Corporation
Convertors/Surgical Specialties
1500 Waukegan Road
McGaw Park, IL 60085
1.800.964.5227
www.cardinal.com/allegiance