



PROTEGRITY

**Triple-layer
protection**

**Advanced
nitrile coating**


**Unparalleled
barrier integrity**

**Precision
grip**

**Effortless
donning**

We've Broken a New Barrier in Barrier Protection.

The Future of Barrier Protection Has Arrived.



It's Protegrity™, the industry's first triple-layer, nitrile-coated, powder-free surgical glove. The Protegrity glove is created with an advanced triple-dip bonding technology that integrates the fit and feel of natural rubber latex with the powerful strength of nitrile for unparalleled barrier integrity.

The natural rubber latex outer layer provides a soft feel, high tensile strength and good elasticity. The middle layer consists of a specially formulated blend of nitrile and natural rubber. This layer enhances the adhesion of the nitrile coating to the latex, enabling the coating to maintain its integrity. The inner layer is a 100% advanced nitrile synthetic coating which provides effortless donning as well as outstanding durability and protection. The combination of these three uniquely bonded layers provides a glove with superior puncture and tear resistance.

Protegrity – it's protection *plus* integrity that equals true peace of mind.

Ordering Information

Catalog Number	Size
2D72N01	5 ½
2D72N11	6
2D72N21	6 ½
2D72N31	7
2D72N41	7 ½
2D72N51	8
2D72N61	8 ½
2D72N71	9

All products packaged
40 pairs per box/5 boxes per case

Professional Services at Allegiance

Our Professional Services Department offers ongoing service and support to ensure customer satisfaction. Our toll-free number provides a direct source for inquiries regarding clinical application and product usage. The Professional Services Group is staffed by registered nurses who offer a full range of clinical support services and product information.

To contact the Professional Services Department, call 1-800-327-7503.

- Unique triple-dipping technology produces bonded layers for optimal barrier protection.
- Inner nitrile synthetic coating facilitates fast and effortless donning.
- Advanced post-forming technology and light chlorination create an exceptional grip, outstanding tactile properties and enhanced comfort.
- Ideal for procedures requiring precision grip.
- Small, protective Allegiance pack is easy to open and saves on storage space.
- Beaded 12" cuff prevents rolldown and provides extended hand and wrist coverage.

To order, call Allegiance Customer Service at 1-800-964-5227.

©Copyright 2001, Allegiance Healthcare Corporation. All rights reserved.
Lit. No. GLV00740 (0201/5M)

♻️ Printed on recycled paper with 10% post-consumer fiber.

 **Allegiance**

a Cardinal Health company

**Allegiance Healthcare Corporation
Gloves**

1500 Waukegan Road
McGaw Park, IL 60085

www.cardinal.com/allegiance