

Medi-Vac®
Tubing and Suction Handles

As the market leader
in clinical suction products,
Medi-Vac offers a
comprehensive line
of sterile suction tubing
and handles to meet
your needs.

 **Allegiance**
a Cardinal Health company

Medi-Vac® Suction Tubing

Medi-Vac sterile suction tubing is a cost-effective solution with the quality and reliability you need to remove fluids from the surgical field.

Features and Benefits

A special formulation provides easy handling and drapability with low coil memory.

Thick tubing walls offer safety and collapse resistance at high vacuum pressures.

Patented Maxi-Grip™ connectors provide a secure fit on a wide range of surgical instruments and equipment.

Available in a variety of formats:

- Convenient precut sterile lengths
- Economical nonsterile 100' rolls
- Standard and custom Medi-Vac canister kits
- Custom Sterile™ kits and PBDS Pathways® kits

The wide variety of sizes and features available from Medi-Vac facilitates economy achieved through standardization while maintaining quality and flexibility.



Bulk Nonsterile Rolls

Tubing is available in 100' rolls with Maxi-Grip connectors at 6' intervals for length customization. A convenient packaging dispenser is provided for easy transport and use.

Cat. Nos. N5100, N6100, N7100



Tube Holding Device

Multipurpose, adhesive-backed holding device is for use with suction tubing or electrical cords. Device allows the tube to slide or secures the tube in a stationary position.

Cat. No. K99



Maxi-Grip™ Connector

Tubing is equipped with Maxi-Grip connectors to accommodate a wide range of surgical instruments and equipment with a secure fit.

Male Connector

In addition to the Maxi-Grip connector, every sterile tube includes a rigid male connector for coupling additional tubes.

Medi-Vac® Suction Handles


Medi-Vac suction handles deliver performance and maneuverability. These lightweight handles are designed to provide optimal control and visibility with minimal tissue trauma.

Features and Benefits

Side holes minimize tissue trauma.

Available in a variety of formats:

- Individually packaged
- Preconnected to patient tubing for easy setup
- Custom Sterile™ kits and PBDS Pathways® kits




Tapered Bulbous Tip

Tapered flanges to reduce tissue trauma with maximum visibility. Available with or without vent.

Cat. Nos. K82, K80

Tip enlarged to show detail.



Bulbous Tip

Rounded to reduce tissue trauma. Available with or without vent.

Cat. Nos. K86V, K86

Tip enlarged to show detail.



Straight Tip

Open tip without unnecessary bulk for better visibility. Available with or without vent.

Cat. Nos. K87V, K87

Tip enlarged to show detail.



Flexi-Clear™ Handle

Provides suction flexibility with clarity and visualization. Available in three tip styles: fine, regular and high. Vented and tubing options available with regular capacity handle.


Cat. Nos. K60, K61, K62, K61V, K61-606, K61-612



Pooler Handle

Multiposition versatility from sump action to pinpoint suctioning. Available with or without tubing.

Cat. Nos. K776, K770



Frazier Handle

Complete with obturator that allows shaping of the shaft to the desired position without occluding. Includes control vent.

Cat. Nos. K71, K72, K73

Medi-Vac® Suction Tubing

- Single Use, Sterile
- Maxi-Grip™ Connectors
- Rigid Male Connector

Cat. No.	Description	Qty. Cs
N52A	5mm (3/16") x 46cm (18")	100
N56A	5mm (3/16") x 1.8m (6')	50
N510	5mm (3/16") x 3.0m (10')	30
N512	5mm (3/16") x 3.7m (12')	30
N520A	5mm (3/16") x 6.1m (20')	25
N66A	6mm (1/4") x 1.8m (6')	45
N610	6mm (1/4") x 3.0m (10')	25
N612	6mm (1/4") x 3.7m (12')	25
N620A	6mm (1/4") x 6.1m (20')	20
N76A	7mm (9/32") x 1.8m (6')	40
N710	7mm (9/32") x 3.0m (10')	20
N712	7mm (9/32") x 3.7m (12')	20
N720A	7mm (9/32") x 6.1m (20')	15

Bulk Suction Tubing

- Single Use, Nonsterile
- Maxi-Grip Connectors Every 6'
- 16 Rigid Male Connectors Included
- Convenient Dispenser Package

Cat. No.	Description	Qty. Cs
N5100	5mm (3/16") x 30.5m (100')	1
N6100	6mm (1/4") x 30.5m (100')	1
N7100	7mm (9/32") x 30.5m (100')	1

Tube Holding Device

- Sterile

Cat. No.	Description	Qty. Cs
K99	Tube/cord holding device	100

All Medi-Vac tubing and suction handle products are 100% latex-free.



Canister Kits Save Time and Money

Medi-Vac preassembled canister kits are the convenient and cost-effective solution for all your suction needs. With canister kits, all items used in the suction process, including canisters, tubing, instrumentation and liquid treatment products, are packaged as a single system. Medi-Vac preassembled kits save time and money by simplifying purchasing, inventory and setup procedures.

For more information or to place an order, contact your Allegiance Healthcare Corporation sales representative, or call 1.800.964.5227.

Yankauer Suction Handles

- Single Use, Sterile

Cat. No.	Description	Qty. Cs
K80	Tapered bulbous tip	50
K82	Tapered bulbous tip, vented	50
K86	Bulbous tip	50
K86V	Bulbous tip, vented	50
K87	Straight tip	50
K87V	Straight tip, vented	50

Yankauer Suction Handles with Preconnected Tubing

- Single Use, Sterile

Cat. No.	Description	Qty. Cs
K83A	Tapered bulbous tip with 6mm (1/4") x 1.8m (6') tubing	50
K84	Tapered bulbous tip with 6mm (1/4") x 2.7m (9') tubing	25
K85	Tapered bulbous tip with 6mm (1/4") x 3.7m (12') tubing	25

Flexi-Clear™ Suction Handles

- Single Use, Sterile

Cat. No.	Description	Qty. Cs
K60	Fine capacity	50
K61	Regular capacity	50
K61V	Regular capacity, vented	50
K61-606	Regular capacity with 6mm (1/4") x 1.8m (6') tubing	50
K61-612	Regular capacity with 6mm (1/4") x 3.7m (12') tubing	25
K62	High capacity	50

Pooler Suction Handles

- Single Use, Sterile
- Sump Action or Pinpoint Suction

Cat. No.	Description	Qty. Cs
K770	Pooler suction handle	50
K776	Pooler suction handle with 6mm (1/4") x 1.8m (6') tubing	20

Frazier Suction Handles

- Single Use, Sterile
- Vented

Cat. No.	Description	Qty. Cs
K71	8 Fr	50
K72	10 Fr	50
K73	12 Fr	50

 **Allegiance**


a Cardinal Health company

**Allegiance Healthcare Corporation
Medi-Vac**

1500 Waukegan Road
McGaw Park, IL 60085

www.cardinal.com/allegiance

©Copyright 2001, Allegiance Healthcare Corporation.
All rights reserved.
Lit. No. MV00450 (0301/5M)

 Printed on recycled paper with 10% post-consumer fiber.