

CREATE YOUR OWN MAXIMUM BARRIER

1 Choose your Altius® RT catheter set.

2 Choose one of the following off-the-shelf insertion packs.



Basic CVC Insertion Pack
(K40174/CVC)



Intermediate CVC Insertion Pack
(K40175/CVC)



Comprehensive Plus CVC
Insertion Pack (K40176/CVC/P1)

OR Create your own insertion pack from our extensive range of components.

ORDERING INFORMATION

ALTIUS® RT ACUTE CATHETERS	
PRODUCT CODE	PRODUCT DESCRIPTION
ART-ADL-1215-K	Altius RT Acute 12Fr 15cm dual lumen haemodialysis catheter
ART-ADL-1220-K	Altius RT Acute 12Fr 20cm dual lumen haemodialysis catheter
ART-ADL-1225-K	Altius RT Acute 12Fr 25cm dual lumen haemodialysis catheter
ART-ADL-1415-K	Altius RT Acute 14Fr 15cm dual lumen haemodialysis catheter
ART-ADL-1420-K	Altius RT Acute 14Fr 20cm dual lumen haemodialysis catheter
ART-ADL-1425-K	Altius RT Acute 14Fr 25cm dual lumen haemodialysis catheter

Find out more at www.kimal.com/AltiusRT

Clinical References

1. Central Venous Dialysis Catheter Dysfunction. Martine Leblanc et al. *Advances in Renal Replacement Therapy*, 1997. 2. A prospective clinical trial to evaluate the microbial barrier of a needleless connector. A.L. Casey et al. *Journal of Hospital Infection* (2007) 65, 212-218. 3. Health Care-Associated Bloodstream Infections Associated with Negative- or Positive-Pressure or Displacement Mechanical Valve Needleless Connectors. William R. Jarvis et al. *Clinical Infectious Diseases* 2009;49:1821-7.

For more information, pricing or samples speak to your local Kimal vascular representative or contact:

UK

Tel: +44 (0)845 437 9542
Fax: +44 (0)845 437 9544
Email: customerservices@kimal.co.uk

Germany

Tel: +49 (0)2133 285 141
Fax: +49 (0)2133 285 1420
Email: info@kimal.de

Middle East

Tel: +97 16 552 8146
Fax: +97 16 552 8147
Email: mea@kimal.co.uk

International

Tel: +44 (0)845 437 9543
Fax: +44 (0)845 437 9544
Email: internationalenquiries@kimal.co.uk

Netherlands

Tel: +49 (0)2133 285 141
Fax: +49 (0)2133 285 1420
Email: info@kimal.nl

www.kimal.com

© 2016 Kimal plc. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal plc. Kimal plc reserve the right to change the specification of its products at any time. Kimal plc, Arundel Road, Uxbridge, UB8 2SA.

L-0203-0011

KIMAL
delivering healthcare innovation

ALTIUS® RT

ACUTE HAEMODIALYSIS

ACUTE RENAL THERAPY

DON'T COMPROMISE. MAXIMISE.

KIMAL
delivering healthcare innovation

MAXIMISE PATIENT SAFETY.


LuerSafe® Connectors

Creates closed system. Reducing the risk of air embolism, blood loss and infections.^{2,3}




360° Rotating Suture Wing

Allows the catheter to rotate after securement.




Multi-Tube™ Technology

Improves quality of the flow rates by reducing turbulence and irregularities.



Intelligent Grooves

Reduce vessel wall occlusions and maximises flow rates.



Soft Flexible PUR

Reduces vessel trauma and eliminates catheter kinking.¹



MAXIMUM PATIENT PROTECTION

- Protection from blood loss ✓
- Protection from air embolism ✓
- Protection from catheter kinking ✓

REDUCED RISK

- Reduces the risk of infection ✓
- Reduces the risk of vessel trauma ✓
- Reduces the risk of vessel occlusions ✓

IMPROVED PERFORMANCE

- Improved catheter patency ✓
- Improved Internal Jugular or Subclavian insertion ✓
- Improved flow rates for critical patients ✓



Stiff, inflexible competitor catheter



Altius RT flexible and durable PUR material