

Medi-Vac® Suction and Fluid Collection Products

Tubing and Suction Handles

An industry leader in fluid management products

Medi-Vac® Suction and Fluid Collection Products offer a comprehensive line of sterile suction tubing and handles to meet your needs.

Our wide variety of sizes and features offer cost-effective fluid management choices to help balance all of your priorities — from safety to cost and quality.



Medi-Vac® Suction Tubing

Medi-Vac® Sterile Suction Tubing is a cost-effective solution with the quality and reliability you need to remove fluids from the surgical field.

- **Special formulation**
Provides easy handling and drapability with low coil memory
- **Thick tubing walls**
Offer safety and collapse-resistance at high vacuum pressures
- **Maxi-Grip™ Connectors**
Patented connectors provide a secure fit on a wide range of surgical instruments and equipment
- **Available in a variety of formats:**
 - Convenient precut sterile lengths
 - Economical, non-sterile 100 ft. rolls
 - Standard and custom Medi-Vac® Canister Kits
 - Presource® Procedure Packs and Presource® PBDS Modules

Bulk non-sterile rolls

Tubing is available in 100 ft. rolls with Maxi-Grip™ Connectors at 6 ft. intervals for length customization. A convenient packaging dispenser is provided for easy transport and use.

Cat. nos. N5100, N6100, N7100

Large-bore suction tubing

Large-bore $\frac{3}{8}$ in. I.D. maximizes flow and drainage capacity. Thick tubing walls offer safety and collapse-resistance.

Cat. no. N1010

Maxi-Grip™ Connector

Tubing is equipped with female Maxi-Grip™ Connectors to accommodate a wide range of surgical instruments and equipment with a secure fit.

Male connector

In addition to the Maxi-Grip™ Connector, every sterile tube includes a rigid male connector for coupling additional tubes.

Canister kits save time and money

Medi-Vac® Preassembled Canister Kits are the convenient and cost-effective solution for all your suction needs. With our canister kits, all items used in the suction process, including canisters, tubing and handles are packaged as a single system. Medi-Vac® Preassembled Kits save time and money by simplifying purchasing, inventory and setup procedures. For more information contact your Cardinal Health sales representative or call 800.964.5227.



Medi-Vac® Suction Handles

Medi-Vac® Suction Handles deliver performance and maneuverability. These lightweight handles are designed to provide optimal control and visibility while helping to minimize tissue trauma.

- **Minimize trauma**
Side holes help minimize tissue trauma
- **Available in a variety of formats:**
 - Individually packaged
 - Preconnected to patient tubing for easy set up
 - Presource® Procedure Packs And Presource® PBDS Module:

Tapered bulbous tip

Tapered flanges to help reduce tissue trauma with maximum visibility. Available with or without vent.

Cat. nos. K82, K80

Bulbous tip

Rounded to help reduce tissue trauma. Available with or without vent.

Cat. nos. K86V, K86

Straight tip

Open tip without unnecessary bulk for good visibility. Available with or without vent.

Cat. nos. K87V, K87

Flexible suction handle

Provides suction flexibility with clarity and visualization. Available in three tip styles: fine, regular and high. Vented and tubing options available with regular-capacity handle.

Cat. nos. K60, K61, K62, K61V, K61-606, K61-612

Poole suction handle


Multiposition versatility from sump action to pinpoint suctioning. Available with or without tubing.

Cat. nos. K776, K770

Frazier suction handle

Complete with obturator that allows shaping of the shaft to the desired position without occluding. Includes control vent.

Cat. nos. K71, K72, K73



Ordering information

All Medi-Vac® Tubing and Suction Handle Products are not made with natural rubber latex.

Tubing and devices

Medi-Vac® Suction Tubing

- Single use, sterile
- Maxi-Grip™ Connectors
- Rigid male connector

Cat. no.	Description	Qty. Cs
N52A	5mm (3/16 in.) x 46cm (18 in.)	100
N56A	5mm (3/16 in.) x 1.8m (6 ft.)	50
N510	5mm (3/16 in.) x 3.0m (10 ft.)	30
N512	5mm (3/16 in.) x 3.7m (12 ft.)	30
N520A	5mm (3/16 in.) x 6.1m (20 ft.)	25
N66A	6mm (1/4 in.) x 1.8m (6 ft.)	45
N610	6mm (1/4 in.) x 3.0m (10 ft.)	25
N612	6mm (1/4 in.) x 3.7m (12 ft.)	25
N620A	6mm (1/4 in.) x 6.1m (20 ft.)	20
N76A	7mm (9/32 in.) x 1.8m (6 ft.)	40
N710	7mm (9/32 in.) x 3.0m (10 ft.)	20
N712	7mm (9/32 in.) x 3.7m (12 ft.)	20
N720A	7mm (9/32 in.) x 6.1m (20 ft.)	15

Large-bore suction tubing

- Single use, sterile
- Large 3/8 in. I.D. maximizes flow and drainage capacity

Cat. no.	Description	Qty. Cs
N1010	10 ft. large-bore tubing	10

Suction handles

Yankauer suction handles

- Single use, sterile

Cat. no.	Description	Qty. Cs
K80	Tapered bulbous tip	50
K82	Tapered bulbous tip, vented	50
K86	Bulbous tip	50
K86V	Bulbous tip, vented	50
K87	Straight tip	50
K87V	Straight tip, vented	50

Yankauer suction handles with preconnected tubing

- Single use, sterile

Cat. no.	Description	Qty. Cs
K83A	Tapered bulbous tip with 6mm (1/4 in.) x 1.8m (6 ft.) tubing	50
K84	Tapered bulbous tip with 6mm (1/4 in.) x 2.7m (9 ft.) tubing	25
K85	Tapered bulbous tip with 6mm (1/4 in.) x 3.7m (12 ft.) tubing	25

Flexible suction handles

- Single use, sterile

Cat. no.	Description	Qty. Cs
K60	Fine capacity	50
K61	Regular capacity	50
K61V	Regular capacity, vented	50
K62	High capacity	50

Pooler suction handles

- Single use, sterile
- Sump action or pinpoint suction

Cat. no.	Description	Qty. Cs
K770	Pooler suction handle	50
K776	Pooler suction handle with 6mm (1/4 in.) x 1.8m (6 ft.) tubing	20

Frazier suction handles

- Single use, sterile
- Vented

Cat. no.	Description	Qty. Cs
K71	8 Fr	50
K72	10 Fr	50
K73	12 Fr	50



Customer Services and Distribution Contact Details

United Kingdom

Tel: +44 (0) 1256365458

Email: uksales@iskushealth.com

Ireland

ROI Tel: +353 1 404 8383

NI Tel: 0845 850 8600

Email: info@iskushealth.com

www.iskushealth.com

CHFM.07.02.2018.1.ISK