

**TECNODRAPE®**

**SURGICAL INCISION  
ANTIMICROBIAL DRAPE**

## **Surgical Infection Prevention Solutions**




**TecnoDrape is effective and essential protection against the contamination of the operating field and surgical incisions. TecnoDrape is designed to prevent the transmission of pathogens from the non-sterile skin towards the sterile surgical field in order to reduce the risk of the patient surgical site infections.**

## SURGICAL INCISION ANTIMICROBIAL DRAPE

**TECNODRAPE** is a surgical incise drape made with **POVIDONE IODINE** and it has a grip tape made of a slow release iodine complex. **TECNODRAPE** is effective protection against the contamination of the operative field and the surgical incision site.

| FEATURES                            | BENEFITS  |
|-------------------------------------|---|
| Iodine Solution                     | Broad spectrum bactericidal action. Fights deep bacterial flora into the top cell layers of the skin.   |
| Gradual release of Iodine particles | Sustained duration of action with the gradual release of iodine molecules for up to six hours from the time of first application.   |
| Acrylate adhesive tape              | Strong adhesion. Perfect adhesion to skin even in very long surgeries. Perfectly isolates the cut borders preventing microorganism's migration. Lock and control of other drapes. No skin damage after removal. Hypoallergenic. Excellent tolerance for the skin. |
| Polyurethane support                | High conformability. Emo-hydro-alcohol repellent. Impenetrable barrier against bacteria. Antiglare. RX Transparent. High extensibility (600%). No memory. Antistatic non-conductor.   |



### PACKAGING

- Single pack with double sterile barrier system
- Easy-opening to avoid cross-contamination
- 2 removable labels for easy traceability
- Internal dispenser for clean storage
- Dimensions and weight are certificated according with Dlgs 81/2008 ex626/94

| NHS SC Code | Product Code | Total Drape Size | Incise Area | Units per Case |
|-------------|--------------|------------------|-------------|----------------|
| VEB961      | IH6635TD     | 15 x 20cm        | 10 x 20cm   | 50             |
| VEB959      | IH6640TD     | 60 x 35cm        | 35 x 35cm   | 50             |
| VEB957      | IH6648TD     | 90 x 60cm        | 60 x 60cm   | 50             |
| VEB960      | IH6650TD     | 90 x 45cm        | 60 x 45cm   | 50             |
| VEB958      | IH6651TD     | 90 x 90cm        | 60 x 90cm   | 24             |

### TECNODRAPE MEDICAL CLASSIFICATION

TECNODRAPE is classified as a Medical Device CE Class III according to directive 93/42/ CEE, notifying body CE0373.

Class III All the devices include a substance that could be considered a medical drug if used separately and by itself, according to article 1 of 65/65/CEE, and that could cause incidental effects to the human body that fall within the III class.



### Customer Services and Distribution Contact Details

#### United Kingdom

Tel: +44 (0) 1256365458

Email: [uksales@iskushealth.com](mailto:uksales@iskushealth.com)

#### Ireland

ROI Tel: +353 1 404 8383

NI Tel: 0845 850 8600

Email: [info@iskushealth.com](mailto:info@iskushealth.com)

[www.iskushealth.com](http://www.iskushealth.com)