

# Compartment Syndrome?

Compass for precise compartment pressure measurement in seconds!



When a crush injury, a car accident, or a sports injury occurs, it raises the question, “Does this patient have compartment syndrome?” Confirming compartment syndrome quickly and confidently can be difficult.

## Compass gives you confidence to make the call

Early detection of compartment syndrome is key to both management and treatment. The sterile, single-use Compass provides a simple, precise, indisputable digital pressure reading of compartment pressure instantly.

### Delays and misdiagnosis can have catastrophic consequences!

Nerves and tissue can die and muscles can be damaged quickly.

With compartment syndrome, early detection is vitally important. For patients who receive an early fasciotomy:

- 16.1% lower rate of amputation
- 7.9% lower rate of infection
- Significantly shorter hospital stays

#### MUSCLE DAMAGE



#### NERVE DAMAGE



**Treatment delays and misdiagnoses are all preventable with the right tool!**



Compass Universal Hg #CUHG1-INT  
Box of 5 Units

- Measures in either mmHg
- Luer lock connection
- Comes with a guidewire port
- Single patient use, sterile, disposable packaging, 5 per box
- Precise to ANSI standards for accurate measurements
- Not made with natural rubber latex

References available upon request.

Compass® is a registered trademark of Mirador Biomedical, Inc. (a subsidiary of Centurion Medical Products Corporation). Centurion® is a registered trademark of Centurion Medical Products Corporation.

©2016 Centurion Medical Products Corporation. All rights reserved. LIT213V2



#### Customer Services and Distribution Contact Details

##### United Kingdom

Tel: +44 (0) 1256365458

Email: [uksales@iskushealth.com](mailto:uksales@iskushealth.com)

##### Ireland

ROI Tel: +353 1 404 8383

NI Tel: 0845 850 8600

Email: [info@iskushealth.com](mailto:info@iskushealth.com)

[www.iskushealth.com](http://www.iskushealth.com)