

## OptiFlo Accessories

### Closed Suction System



#### Versatile system

The OptiFlo closed suction system is equipped for almost all indications. OptiFlo is available in a great number of sizes and lengths, with various durations of service life and special features.

#### Tried and tested solutions

You also have access to our range of accessories. This is continuously being enhanced by new products derived from knowledge gained from our experience in the field to assist you in everyday clinic work.

#### Safe partner

We are a partner that you can depend upon. Do have any questions or do you require a further feature? Just get in touch with us.



Closed patient suction system



Endotracheal suction tube in use

# OptiFlo Accessories

One system for many indications



Catheter mount, smooth bore tube

## Catheter mount, smooth bore tube, 15M–15F

Length	REF	Packing unit
5 cm	13.502.00.005	25
10 cm	43.500.15.100	25
15 cm	13.502.00.015	25



Catheter mount, extendable

## Catheter mount, extendable, 15M–15F

Length	REF	Packing unit
5 cm	13.505.00.005	25
10 cm	13.505.22.010	30
15 cm	13.505.00.015	25



Tracheal suction system

## Tracheal suction system

Description	REF	Packing unit
Female connector / fingertip	07.098.00.000	180
Male connector / female connector	07.098.00.003	180



Endotracheal suction tube

## Endotracheal suction tube

Description	REF	Packing unit
Vacuum control / 2 clamps	05.000.22.657	30
Vacuum control / 2 clamps / direct OptiFlo connection	05.000.22.658	30

## Rinsing solution

Filling capacity	REF	Packing unit
15 ml	43.004.70.110	144

## Blue adapter for OptiFlo

	REF	Packing unit
	43.004.60.099	10

## Customer Services and Distribution Contact Details

### United Kingdom

Tel: +44 (0) 1256365458

Email: [uksales@iskushealth.com](mailto:uksales@iskushealth.com)

### Ireland

ROI Tel: +353 1 404 8383

NI Tel: 0845 850 8600

Email: [info@iskushealth.com](mailto:info@iskushealth.com)

[www.iskushealth.com](http://www.iskushealth.com)