



ISKUSHEALTH


IMPROVING PATIENT OUTCOMES

Needleless Connectors

DELTA P. VALVE



**The smallest and lightest neutral
needle free connector with the
highest flow rate**



Neutral needle free connector with dedicated fluid path

The “DELTA P. VALVE”

- Priming volume: 0,04 ml
- High flow rates: 160 ml/min
- Light weight: 1,5 grams - suitable for neonatal patients
- Neutral displacement for catheter tip patency
- Split septum system without any internal mechanical part
- Totally transparent: ability to visualise internal housing for complete flushing of the device
- Easy to disinfect the surface of the device



Iskus Health UK Ltd

Unit 5, Swallowfield Way, Hayes,
Middlesex, London, UB3 1DQ,
United Kingdom

Phone : +44 (0) 208 869 6509

E-mail: uksales@iskushealth.com

www.iskushealth.com

Manufactured by:

 **DELTA MED**
MEDICAL DEVICES
Solutions for Health

DELTA MED SPA, Viadana (MN) - Italy