


# FlexAgris™ Catheter Securement Devices



Designed for


Securing Catheters with  
Comfort



## FlexAgrip™ Catheter Securement Device

The only catheter securement device with shape memory technology

The patented FlexAgrip™ Catheter Securement Device, with shape memory properties, was designed to take the shape of any type of catheter providing superior catheter securement. The FlexAgrip™ device was engineered to withstand high pull forces while providing patient comfort due to its soft, breathable materials used.



➔ **Increased catheter dwell time**

➔ **Reduced catheter complications**

➔ **Less application steps**

➔ **Compliments existing IV dressing protocols**


Product Code: FGS0109  
NHS SC: ELW854



## Pharmapore<sup>®</sup> Catheter Securement Dressing

The only securement dressing with shape memory technology





The patented all-in-one Pharmapore<sup>®</sup> Catheter Securement Dressing provides insertion site protection while securing the catheter. With FlexAgrip being integrated within the dressing, it allows the dressing to take the shape of any type of catheter providing superior catheter securement while providing patient comfort due to its soft, breathable materials used.



- All-in-one catheter securement dressing providing insertion site protection & stabilizing the catheter
- Large viewing window allows for constant visual inspection of insertion site
- Breathable, waterproof barrier protecting the skin from outside contaminants
- Designed to be worn for up to 7 days
- Sutureless securement

Product Code: CSD0422  
NHS SC: ELY855

**FlexAgrip™** catheter securement devices & dressings come in a variety of types and sizes that support a wide range of IV applications.

Product	Product Code	Product Description	Unique Features & Catheter Application	Size (cm)	Dressings/Box (Boxes/Case)
	CSD0422 NHS SC ELY855	Pharmapore® Catheter Securement Dressing	FlexAgrip™ integrated into an occlusive dressing. Applied to CVC/PICC	9.75 x 11.75 4.0 x 9.0	50 (24)
	CSD0119 NHS SC ELY853	Pharmapore® Catheter Securement Dressing	FlexAgrip™ integrated into an occlusive dressing plus securement device for additional securement. Applied to CVC/PICC	9.75 x 11.75 4.0 x 9.0 4.0 x 8.5	50 (24)
	CSD1129 NHS SC ELW1003	Pharmapore® Catheter Securement Dressing	FlexAgrip™ integrated into an occlusive dressing. Applied to PIV w/ Rotating Hubs	7.0 x 9.0 2.25 x 6.6	50 (24)
	FGS0109 NHSSC ELW854	FlexAgrip Catheter Securement Device	FlexAgrip™ technology applied to PIV or arterial catheters	4.0 x 8.5	100 (16)



**Tel:** +44 (0) 208 869 6509

**Email:** [uksales@iskushealth.com](mailto:uksales@iskushealth.com)

[www.iskushealth.com](http://www.iskushealth.com)

FlexAgrip™ & Pharmapore® Catheter Securement Devices