

PROTECTION

MARLUX with Microban SilverShield helps to keep staff and patients safer



Pathogens can survive on surfaces for long periods of time and are usually invisible.



Cross-contamination between staff, patients, family and surfaces could increase the risk of hospital acquired infections (HAI's).



Microban SilverShield is proven to be bactericidal and kills up to 99.99% of bacteria.¹



Microban SilverShield is built-in at the point of manufacture to ensure an even distribution.

Untreated polypropylene

May inhibit growth but does not kill bacteria

risk of transfer

Coated polypropylene

Cannot ensure even distribution and could wear off or wash off


risk of transfer

Traditional curtains

Invisible pathogens could build up between laundry rounds

risk of transfer


Unprotected



2 hrs

24 hrs

with Microban



2 hrs

24 hrs



SAFETY

Fire Retardant

MARLUX curtains are manufactured and stocked in the UK and follow a rigorous batch testing regime.




Our Fire Retardant testing is carried out in the UK by an accredited company for every 5,000m of curtain produced and every curtain is batch coded for traceability.

Fire retardant to BS 5867-2: 2008 type C

Manufacturing


MARLUX has a manufacturing facility in Coleshill (UK), where our production team take large rolls of raw material through cutting, pleating, and ultrasonic welding to create a finished product ready for worldwide distribution.

MARLUX curtains are made from 100% polypropylene which is 100% recyclable.



CONVENIENCE

Fittings for any track



FastFit



Resistoid



Silent Gliss + Mark 4 insert




Goelst + Mark 4 insert



Mova Track + Mark 5 insert



Antiference + Mark 5 insert



Uni-Glide




Silent Gliss



Goelst



Antiference



Universal



Resistoid + FFA-hook

Universal can fit onto existing hooks already in place.

Alternatively you can use the hooks below to convert them.



Mark 4



Mark 5



FFA-hook



UVS-hook



UGN-hook

PRODUCT RANGE

ORDER CODES



MARLUX curtains are available in 1.95m height and 100gsm.



FastFit



Universal



Uni-Glide



Printed

Available in the following colours:



Pastel Blue



Bali Sand



Forest Green



Summer Blue

Width	Colour	Box qty	FastFit	Universal	Uni-Glide
7.2m	Summer Blue	8	FSB72/1950/A	USB72/1950/A	GSB72/1950/A
7.2m	Forest Green	8	FFG72/1950/A	UFG72/1950/A	GFG72/1950/A
7.2m	Bali Sand	8	FBS72/1950/A	UBS72/1950/A	GBS72/1950/A
7.2m	Pastel Blue	8	FPB72/1950/A	UPB72/1950/A	GPB72/1950/A
4.2m	Summer Blue	12	FSB42/1950/A	USB42/1950/A	GSB42/1950/A
4.2m	Forest Green	12	FFG42/1950/A	UFG42/1950/A	GFG42/1950/A
4.2m	Bali Sand	12	FBS42/1950/A	UBS42/1950/A	GBS42/1950/A
4.2m	Pastel Blue	12	FPB42/1950/A	UPB42/1950/A	GPB42/1950/A
1.8m	Summer Blue	32	FSB18/1950/A	USB18/1950/A	GSB18/1950/A
1.8m	Forest Green	32	FFG18/1950/A	UFG18/1950/A	GFG18/1950/A
1.8m	Bali Sand	32	FBS18/1950/A	UBS18/1950/A	GBS18/1950/A
1.8m	Pastel Blue	32	FPB18/1950/A	UPB18/1950/A	GPB18/1950/A

Product	Code	Box qty
Universal to Fast Fit hooks	FFA-hook	200
Universal to Uni-Glide lollipop hooks	UGN-hook	200
Universal Seahorse Hook	UVS-hook	200
Adapta-Track MK4 insert	MK4	18m/box
Adapta-Track MK5 insert	MK5	18m/box
Silicone wipes	SW/100	100

When to change
 The Care Quality Commission guidance stipulates that disposable privacy curtains should be changed at least once every 6 months. The guidance from both the Infection Control Nurses Association and the National Patient Safety Agency reiterate the same responsibility. Soiled curtains should be changed immediately.

Custom sizes or printed curtains available on request
REFERENCES
 *Industrial microbiological services Ltd. IB10-5083-100 October 2012

MARLUX MEDICAL
 MARLUX MEDICAL Ltd
 Business Unit
 Part of OrthoD Group Ltd
 Industrial Park
 Bourton-on-the-Water
 Gloucestershire
 GL54 2HQ
 United Kingdom

tel. +44 (0) 121 783 5777
 email. sales@marlux.co.uk
 web. www.marlux.co.uk



MARLUX DISPOSABLE PRIVACY CURTAINS

