

Woman Centered Caesarean Birth Drape



**Delivering immediate
"Skin to Skin" contact
between baby & mother!**

A Woman Centered C-Section Drape

An aperture in the drape allows the baby to pass through following its delivery for immediate skin to skin contact with the mother, which stimulates bonding and breast feeding^{1,2}.

Encourages and facilitates the inclusion of the father/partner during the birth for the cutting of the cord^{1,2}.

Overall enhances the experience for mothers creating a more Gentle Caesarean Section birth.

Parents feedback²


"to be the first to hold my baby was so special..."

"to be included as a Dad and get to cut the cord was very important to me..."

"I felt so relaxed and calm, a lovely experience..."



This new innovative Woman Centered Caesarean Drape was developed in collaboration with obstetricians and midwives to support the current research that early "skin to skin" contact between a new born baby and mother is hugely beneficial in enhancing the birthing experience for those women which require a C-Section delivery. #ImprovingPatientOutcomes



Features of the drape include:

30 x 35cm Baby Aperture with closeable flap

360° fluid control pouch with contoured foam ring, Cord-Holding Tabs and double drainage ports

Armboard covers

30cm round fenestration with clear incise

Latex free and ETO sterile

Ordering information

Product code	NHS PC	Description	Case Qty
MTT0523-ISK	VJL280	Cesarean drape with baby transfer aperture (sterile)	8

References

- Hunter et al. BMC Pregnancy and Childbirth (2017) 17:322
- O'Sullivan et al. Poster; University Maternity Hospital Limerick, Woman Centred Approach to Caesarean Birth (2017)



Customer Services and Distribution Contact Details

United Kingdom

Tel: +44 (0) 1256 365458

Email: uksales@iskushealth.com

Ireland

ROI Tel: +353 1 404 8383

NI Tel: 0845 850 8600

Email: info@iskushealth.com

www.iskushealth.com