



# FFP3 RESPIRATOR FACE MASKS



Give yourself peace of mind: adapted for all your needs of protection against airbourn particles.

## SPECIFICATIONS

- » Non-valved flat fold
- » Single use
- » Individually packed for improved hygiene
- » Concave shape with horizontal fold *for a better seal*
- » Comfortable fit and *large dimension*
- » Made of non-woven
- » 2 soft elastic head straps
- » 1 nose clip
- » Filtration and breathability insured  
*Filtration ≥ 99% (FFP3) and very low breathing resistance; As low as maximums for FFP1*
- » Made and certified in France
- » Stock held in UK

## LOGISTICS

	Packaging	Reference
FFP3	Case of 4 boxes of 25 masks	1330-C
	Case of 200 masks in bulk	1330-V



## REGULATION

- » Individual Protective Equipments cat III compliant to PPE Regulation 2016/425
- » Filtering half masks to protect against particles compliant to EN 149: 2001+ A1: 2009FFP3: Small particles filtration ≥ 99%



UK Tel: +44 (0) 208 869 6509  
 Email: [uksales@iskushealth.com](mailto:uksales@iskushealth.com)  
[www.iskushealth.com](http://www.iskushealth.com)

