

Surgical Drapes, Gowns and Ancillary Products



CardinalHealth
Essential to care™


Distributed by



ISKUS HEALTH
IMPROVING PATIENT OUTCOMES

Cardinal Health provides a complete selection of surgical drapes, gowns and ancillary products that offers patients and clinicians the protection they want and need.

Our enhanced line of Tiburon® Surgical Drapes provides improved drapeability*, additional cord management features* and AAMI level 4 protection. Additionally, our surgical gowns offer AAMI level 3 (Astound® Surgical Gowns) and AAMI level 4 protection (SmartGown™ Surgical Gowns) to help keep you safe from potentially harmful elements, allowing you to focus on what matters most – your patients. In developing our portfolio of products, we worked closely with healthcare professionals, like you, to make sure that your needs and expectations were not just met, but exceeded.



*As compared to prior Tiburon® Surgical Drape design.

AAMI Levels determined per ANSI/AAMI PB70:2012

Contents

Surgical gowns.....	5
SmartGown™ Surgical Gowns.....	5
Astound® Surgical Gowns.....	5
Sterile sleeve	5
Tiburon® Surgical Drapes and Packs.....	6-59
Ancillary.....	7-12
Ancillary patient drapes.....	8-10
Drape accessories	11
Equipment drapes and covers.....	12
Angiography	13-15
Cardiovascular	16-18
Central Line.....	19-20
Cystoscopy	21-22
General surgery.....	23-26
Neurosurgery.....	27-28
Obstetrics	29-31
Ophthalmology and ENT	32-37
Orthopedics.....	38-44
Pacemaker	45-46
Pediatrics	47-48
Peri-GYN	49-50
Special procedures	51-52
Standard packs.....	53-59
Indexes.....	60-65
Sterile number index.....	60-61
Nonsterile number index.....	62-63
Product index.....	64-65



Surgical gowns

Be in the know and know you are protected with Cardinal Health™ Surgical Gowns.

Inspired by long days on your feet, our surgical gowns stand up to the rigors of your work, helping to keep you safe from potentially harmful elements.

In creating SmartGown™ and Astound® Surgical Gowns, we worked closely with healthcare professionals like you to make sure your needs and expectations were not just met, but exceeded.



What will you be wearing?

It's important to choose the right gown for the right procedure. While there are many factors that go into choosing the right surgical gown, the questions below may help as you evaluate the different gowns and levels of protection available:

- What type of procedure is being preformed?
- How long will the procedure last?
- How much fluid will be present during the case?
- What is the clinician's role in the procedure?²



All Cardinal Health™ Surgical Drapes and Gowns meet the EN13795 Standard

AAMI Levels determined per ANSI/AAMI PB70:2012

SmartGown™ Surgical Gowns



SmartGown™ Surgical Gowns. Fully Impervious

Sterile No.	Type	Qty. Cs
89015	Large	Set-in Sleeve 20
89045	X-Large	Set-in Sleeve 18
89049	X-Large X-Long	Set-in Sleeve 14
39049	X-Large, X-Long	Raglan Sleeve 14
39075	XX-Large	Raglan Sleeve 14
32474	X-Large, X-Long, A-line	Raglan Sleeve 16
39079	XX-Large X-Long	Raglan Sleeve 14

Surgical Gowns

89013 SmartGown™ Surgical Gowns
and 890015N SmartGown™ Surgical
Gowns and Scrub Nurse

Astound® Surgical Gowns



Astound® Standard Gowns

Sterile No.	Size	Type	Qty. Cs
ASG9515	Large	Set-in Sleeve	20
ASG9545	X-Large	Set-in Sleeve	20
ASG9575	XX-Large	Set-in Sleeve	18


Astound® Fabric-Reinforced Gowns

Sterile No.	Size	Type	Qty. Cs
ASG9511	Large	Set-in Sleeve	20
ASG9541	X-Large	Set-in Sleeve	18
ASG9541EL	X-Large, X-Long	Set-in Sleeve	18

Astound® Poly-Reinforced Gowns

Sterile No.	Size	Type	Qty. Cs
9010	Large	Set-in Sleeve	20
9040	X-Large	Set-in Sleeve	20
9070	XX-Large	Set-in Sleeve	18
9041EL	X-Large X Long	Set-in Sleeve	18

Sterile Sleeve



Sterile Disposable Sleeve

- Quick sterile coverage for contaminated gown sleeves
- Elastic at both ends for ease of donning
- 21.5 in. minimum length

Sterile No.	Qty. Bx	Qty. Cs
599	30	60/cs

Tiburon® Surgical Drapes and Packs

Tiburon® Surgical Drapes are engineered to pass the industry standard tests that measure fluid penetration to provide AAMI level 4 protection. They are resistant to abrasion and puncture, yet their flexible, cloth-like properties allow them to drape naturally, following the contours of the patient.

Classified as a Level 4 device under the AAMI guidelines for barrier performance in the critical zone, Tiburon® Surgical Drapes provide the highest level of barrier protection available in a surgical drape.

You asked. We Listened.

The next-generation of Tiburon® Surgical Drapes delivers improved drapeability and additional cord management features.

Features and benefits of the next-generation Tiburon® Surgical Drape include:*

- AAMI Level 4 protection
- High-performance SMS material that customers know and trust
- Designed to stay in place throughout the procedure; avoids slippage
- Softer fabric for enhanced drapeability
- Velcro® Tabs for cord management on select codes
- Fabric anesthesia screen on select codes




Enhanced by technology

- We're pleased to bring to the market surgical drapes with proven SMS fabric that continue to deliver protection and performance.



Improved by you

- You asked for a more drapeable fabric; we listened.
- You asked for better cord management; we delivered.
- Your feedback helped us develop the next-generation Tiburon® Surgical Drapes.



Tiburon® Surgical Drapes and Packs


Ancillary

Ancillary Patient Drapes



Ancillary

18 in. (45cm)




D1000

24 in. (61cm)



D1089

44 in. (112cm)




9462 Tiburon™ Utility Drape Sheet

Drape Towels

Sterile No.	Nonsterile No.	Description	Qty. Cs
D1000	D1000NS	Towel drape 18 x 12 in. (45 x 30cm), Adhesive 18 x 1.5 in. (45 x 3.8cm)	10
D1010	D1010NS	Towel drape 18 x 24 in. (45 x 61cm), Adhesive 18 x 1.5 in. (45 x 3.8cm)	10
D1089	D1089NS	Utility drape (2 drapes per pkg) 24 x 26 in. (61 x 66cm)	20 pkgs/cs

Drape Sheets

Sterile No.	Nonsterile No.	Description	Qty. Cs
9358	9358N	Tiburon® Half Sheet, 40 x 57 in. (101 x 145cm)	20
9355	9355N	Tiburon® Medium Sheet, 40 x 71 in. (101 x 180cm)	20
9349	9349N	Tiburon® Three-Quarter Sheet, 52 x 76 in. (132 x 193cm)	20
9446	9446N	Tiburon® Large Sheet, 72 x 85 in. (183 x 216cm)	30
9444	9444N	Tiburon® X-Large Sheet, 77 x 98 in. (195 x 249cm)	12
9462	9462N	Tiburon® Utility Drape Sheet, 44 x 60 in. (112 x 152cm), 4 x 6 in. (10 x 15cm) poly-reinforced fenestration	20



Tiburon® Surgical Drapes and Packs
Ancillary Patient Drapes




Ancillary

Drape Sheet, Square-Folded


Sterile No.	Nonsterile No.	Description	Qty. Cs
9461	9461N	Tiburon® Drape Sheet, square-folded, 38.5 x 38.5 in. (98 x 98cm)	30

Under-Buttocks Drapes


Sterile No.	Nonsterile No.	Description	Qty. Cs
8482	8482N	Under-buttocks drape, with fluid control pouch, filtering screen, drainage port and reinforcement	20
8483	8483N	Under-buttocks drape, reinforced poly, 9 in. cuff, 36 x 49 in. (91 x 125cm)	40
8484	8484N	Under-buttocks drape with fluid control pouch, reinforcement	20
8487	8487N	Under-buttocks drape with graduations, fluid control pouch with filtering screen and drainage port, plastic backing	20



Ancillary Patient Drapes




Ancillary



8430

Leggings

Sterile No.	Nonsterile No.	Description	Qty. Cs
8420	8420NA	Leggings, folded, 6 in. cuff, 28 x 43 in. (71 x 109cm), 2/pkg.	20
8430	8430NA	Clear leggings, folded, 6 in. cuff, 31 x 55 in. (79 x 140cm), 2/pkg	20



D1015

Poly U-Drapes

Sterile No.	Nonsterile No.	Description	Size	Qty. Cs
8475	8475N	Blue poly U-drape, 8 x 30 in. split with adhesive	60 x 84 in. (152 x 213cm)	22
8476	8476NA	Blue poly U-drape, 6 x 21 in. split with adhesive	60 x 72 in. (152 x 183cm)	24
D1015	D1015NS	Clear poly U-drape, Adhesive 9.5 x 19.5 in. (24 x 50cm), split 6 x 19 in. (15 x 48cm)	48 x 52 in. (121 x 132cm)	10
D1019	D1019NS	Clear poly U-drape, Adhesive 9.5 x 19.5 in. (24 x 50cm), Adhesive tab 3 x 1 in. (7 x 2.5cm), split 6 x 19 in. (15 x 48cm)	48 x 71 in. (121 x 180cm)	10




Tiburon® Surgical Drapes and Packs Drape Accessories

Ancillary

Adhesives and Drape Tape

Sterile No.	Nonsterile No.	Description	Qty. Cs
3524	N/A	Drape tape, 4 x 18.5 in. (10 x 47cm)	2

4 in. (10cm)




3524

Pouches

Sterile No.	Nonsterile No.	Description	Qty. Cs
D1016	D1016NS	Instrument pouch 6.75 x 11.5 in. (17 x 29cm), Adhesive (2) 6.75 x 2 in. (17 x 5cm)	10
D1018L	D1018LNS	Instrument pouch 10 x 18 in. (26 x 45cm), Adhesive (2) 10 x 2 in. (26 x 5cm)	10

20 in. (51cm)




D1016

Pouches

Sterile No.	Nonsterile No.	Description	Qty. Cs
D1018	D1018NS	Fluid pouch 20 x 19 in. (51 x 48cm), Adhesive 20 x 3.75 in. (50 x 9cm)	10
D1072	D1072NS	Aperture pouch 18 x 9 in. (45 x 22cm), Aperture 4.75 x 3.75 in. (12 x 9cm), Adhesive (9.5 x 9 in.)	10

*Denotes a 6-8 week lead time

6.75 in. (17cm)



D1018

Isolation Bags

Sterile No.	Nonsterile No.	Description	Qty. Cs
D1003	D1003NS	Isolation bag 20 x 20 in. (51 x 51cm)	10

Equipment Drapes and Covers

Standard Back Tables Covers


Sterile No.	Nonsterile No.	Description	Qty. Cs
8376	8145N	Standard back table cover, 44 x 75 in. (112 x 190cm), SMS reinforcement, 23 x 75 in. (58 x 190cm)	22
8377	8146N	Standard back table cover, 44 x 90 in. (112 x 228cm), SMS reinforcement, 23 x 90 in. (58 x 228cm)	22
8378	8184N	Standard back table cover, 50 x 90 in. (127 x 228cm), SMS reinforcement, 29 x 90 in. (73 x 228cm)	24

Heavy-Duty Back Table Covers

Sterile No.	Nonsterile No.	Description	Qty. Cs
8186	8186N	Heavy-duty back table cover, 60 x 90 in. (152 x 229cm), thick SMS reinforcement, 29 x 90 in. (74 x 229cm)	24

Mayo Stand Covers


Sterile No.	Nonsterile No.	Description	Qty. Cs
8337	8337N	Mayo stand cover, reinforced poly, 23 x 55 in. (58 x 140cm)	30
8339	8339NI	Mayo stand cover, reinforced poly, 29.5 x 55 in. (75 x 140cm)	22




D1012

Equipment Covers

Sterile No.	Nonsterile No.	Description	Qty. Cs
D1012	D1012NS	Fluoroscope Cover 36 x 44 in. (91 x 111cm), Adhesive strips 36 x 0.5 in. (91 x 1.2cm), Adhesive tabs 2.5 x 4.25 in. (6 x 11cm)	10
D2158	D2158NS	Instrument Cover 24 x 18 in. (60 x 45cm), Adhesive 24 x 0.5 in. (60 x 1.2cm)	10
D2008A	D2008NS	Clear Cover 44 x 72 in. (111 x 182cm)	25



D2158



Tiburon® Surgical Drapes and Packs

Angiography

Tiburon® Surgical Drapes and Packs


Angiography



Angiography

Tiburon® Femoral Angiography Drape


Sterile No.	Nonsterile No.	Qty. Cs
29460	29460N	16



- Absorbent reinforcement, 28.75 x 44 in. (73 x 112cm)
- Catheter pouch, 8.5 x 11.5 in. (21 x 29cm)
- Clear side panels, 27.5 x 112 in. (70 x 285cm)
- Fenestration, dual circular, 6 in. (15cm), with circular incise aperture, 3.5 in. (9cm)

Extra Large Femoral Angiography Drape


Sterile No.	Nonsterile No.	Qty. Cs
29570	29570N	8



- Absorbent reinforcement, 33.5 x 135 in. (85 x 343cm)
- Clear foot panel, 80 x 25 in. (203 x 63.5cm)
- Clear side panels, 23.25 x 135 in. (60 x 343cm)
- Fenestration, dual circular, 6 in. (15cm), with circular incise aperture, 3.5 in. (9cm)

Tiburon® Femoral Angiography Drape


Sterile No.	Nonsterile No.	Qty. Cs
29456	29456N	14



- Absorbent reinforcement, 29 x 44 in. (74 x 112cm)
- Fenestration, dual circular, 6 in. (15cm), with circular incise aperture, 3.5 in. (9cm)
- Single clear side panel on patient's right side, 26.5 x 112 in. (67 x 284cm)

Angiography Sheet


Sterile No.	Nonsterile No.	Qty. Cs
D2244	D2244NS	20



- Apertures: 3.5 x 4.5 in. (9 x 11cm)

Angiography Sheet


Sterile No.	Nonsterile No.	Qty. Cs
D1195	D1195NS	20




- Apertures: 3.5 x 4.5 in. (9 x 11cm)
- Moldable strips on right panel

Radial/Femoral Angiography Drape

Sterile No.	Nonsterile No.	Qty. Cs
29529	29529N	12



- Absorbent reinforcement
- Clear side and foot panels
- Fenestration, dual circular, 6 in. (15cm), with circular incise aperture, 3.5 in. (9cm)
- Malleable bars
- Radial access on patient right, 2 x 3 in. (5 x 7.5cm) oval incise aperture



Tiburon® Surgical Drapes and Packs

Angiography

Tiburon® Brachial Angiography Drape

Sterile No.	Nonsterile No.	Qty. Cs
29458	29458N	40



38 in. (96cm)

44 in. (112cm)

- Absorbent reinforcement,
12 x 12 in. (30 x 30cm)
- Fenestration, circular incise aperture,
4 in. (10cm)

Angiography



Tiburon® Surgical Drapes and Packs


Cardiovascular

Tiburon® Surgical Drapes and Packs
Cardiovascular

Cardiovascular

Tiburon® Cardiovascular Split Drape with Pouches


Sterile No.	Nonsterile No.	Qty. Cs
9155	9155N	12



- Absorbent reinforcement, 10 x 57.5 in. (25.4 x 146cm) (sides) and 20 x 15 in. (51 x 38cm) (top)
- Armboard covers
- Clear anesthesia screen, 100 x 28.5 in. (254 x 72cm)
- Clear pouches on patient right and left
- Clipping tab
- Split with adhesive, 7 x 81 in. (18 x 206cm)

Tiburon® Cardiovascular Incise Drape


Sterile No.	Nonsterile No.	Qty. Cs
9156	9156N	7



- Non-absorbent reinforcement
- Armboard covers
- Clipping tabs
- Fenestration, chest incise, 18 x 22 in. (46 x 56cm)
- Fenestration, leg incise, 28 x 45 in. (71 x 114cm)

Tiburon® Cardiovascular Split Drape


Sterile No.	Nonsterile No.	Qty. Cs
29154	29154N	8



- Absorbent reinforcement
- Armboard covers
- Split with adhesive, 5 x 80 in. (13 x 203cm)

Tiburon® Cardiovascular Peri-Groin Drape


Sterile No.	Nonsterile No.	Qty. Cs
9154	9154N	10



- Absorbent reinforcement, 30 x 42 in. (76 x 106cm)
- Attached groin sheet, 36 x 11 in. (91 x 28cm)
- Splits, dual with adhesive, 6 x 20 in. (15 x 51cm)
- Included in pack 29158

Tiburon® Cardiovascular Incise Drape

Sterile No.	Nonsterile No.	Qty. Cs
9162	9162N	7



- Absorbent reinforcement, 42 x 88 in. (107 x 223cm)
- Armboard covers
- Fenestration, chest incise, 14 x 18 in. (25.5 x 46cm)
- Fenestration, leg incise, 24 x 41 in. (61 x 104cm)
- Instrument troughs, 8 x 88 in. (20 x 223cm)
- Velcro® Tabs


Tiburon® Surgical Drapes and Packs

Cardiovascular




Cardiovascular

Tiburon® Head, Side and Foot Drapes

Sterile No.	Nonsterile No.	Qty. Cs
N/A	29462N	
N/A	29467N	
N/A	29468N 99 in. (251cm)	
		
		
		

- Absorbent reinforcement
- Cord-holding tabs
- Included in pack 29153

Tiburon® Bottom Drape with Reinforcement

Sterile No.	Nonsterile No.	Qty. Cs
N/A	29464N	
		

- Absorbent reinforcement, 26 x 30 in. (66 x 76cm)

Cardiovascular Packs

Nonsterile No.	Description	29153 Cardiovascular Pack with reinforced Gowns
9150N	Tiburon® Cardiovascular Drape	
29467N	Tiburon® Head Drape, 99 x 53 in. (251 x 134.6cm), adhesive, absorbent reinforcement	1
29462N	Tiburon® Side Drape, 38 x 77 in. (96 x 195cm), adhesive, absorbent reinforcement	2
29468N	Tiburon® Foot Sheet, 75 x 53 in. (190 x 134.6cm), adhesive, absorbent reinforcement	1
9349N	Tiburon® Three-Quarter Drape Sheet, 52 x 76 in. (132 x 193cm)	
9154N	Tiburon® Cardiovascular Peri-Groin Drape	
9158N	Tiburon® Cardiovascular Split Drape II	
29464N	Tiburon® Bottom Drape with reinforcement	1
9541N	Astound® Gown, sterile back, fabric reinforced, X-large	
9040N	Astound® Gown, poly reinforced, X-large	2
90780N	Astound® Gown, reinforced, X-large, with towel in outer wrap	1
90680NA	Astound® Gown, reinforced, large, with towel in outer wrap	
8520N	Absorbent towel	1
7554N	Tiburon® Drape Towel, with adhesive strip, 19 x 30 in. (48 x 76cm)	
8694-60	Suture bag	1
8337N	Mayo stand cover, reinforced poly, 23 x 55 in. (58 x 140cm)	
8188N	Overhead table cover, 80 x 90 in. (203 x 228.6cm), reinforced poly with a 44 x 72 in. (112 x 183cm) fabric sterile field extension, SMS reinforcement, 40 x 90 in. (101.6 x 228.6cm)	1
8554N	Non-absorbent drape towel with adhesive strip, 15 x 26 in. (38 x 66cm)	6
9460N	Peri-towel	1
8470N	Plastic sheet, 50 x 59 in. (127 x 150cm)	1
8152N	Heavy-duty back table cover, 80 x 110 in. (203 x 279cm), SMS reinforcement, 37 x 110 in. (94 x 279cm)	
8143N	Standard back table cover, 50 x 90 in. (127 x 228.6cm), SMS reinforcement, 29 x 90 in. (73.6 x 228.6cm)	1
	Quantity per case	3





Tiburon® Surgical Drapes and Packs

Central line


Central line

Central line

Full Body Central Line Insertion Drape

Sterile No.	Nonsterile No.	Qty. Cs	Sterile No.	Nonsterile No.	Qty. Cs
29322	29322N	10	29527	29527N	14
 <p>77 in. (195cm) 108 in. (275.5cm)</p>					
 <p>77.5 in. (197cm) 120 in. (305cm)</p>					

- Clear window, 18 x 18 in. (46 x 46cm)
- Fenestration with adhesive, 6 x 6 in. (15 x 15cm)
- Synthetic absorbent reinforcement, 41.5 x 41.5 in. (105 x 105cm)
- Absorbent reinforcement, 36 x 57.5 in. (91 x 146cm)
- Fenestrations with incise aperture allow for intrajugular or subclavian access
- Clear side panels, 20.75 x 120 in. (53 x 305cm)
- Easy removal feature on both fenestrations



Tiburon® Surgical Drapes and Packs

Cystoscopy

Tiburon® Surgical Drapes and Packs


Cystoscopy



Cystoscopy

Tiburon® Cystoscopy "T" Drape


Sterile No. N/A **Nonsterile No.** 29437NA **Qty. Cs**



- Absorbent reinforcement, 12 x 21 in. (30.5 x 53cm)
- Cord-holding tabs
- Fenestration, 3 x 6 in. (7.5 x 15cm)
- Mesh screen, 10 x 10 in. (25 x 25cm)
- Included in packs 29470, 29471, 29480

Tiburon® Urological Drape


Sterile No. N/A **Nonsterile No.** 29451N **Qty. Cs**



- Absorbent reinforcement, 12 x 19 in. (30.5 x 48cm)
- Attached leggings
- Cord-holding tabs
- Fenestration, 3 x 6 in. (7.5 x 15cm)
- Mesh screen, 14 x 8 in. (35.5 x 20cm)
- Included in pack 29360


Tiburon® Lithotomy Drape with Fluid Control Pouch

Sterile No. 29455 **Nonsterile No.** 29455NA **Qty. Cs** 8



- Absorbent reinforcement, 20 x 19 in. (51 x 48cm)
- Anesthesia screen
- Attached leggings
- Cord-holding tabs
- Fenestration with adhesive, 2.5 x 6.5 in. (6 x 16.5cm)
- Fluid control pouch with filtering screen, suction port, drawstring feature and side trough
- Mesh screen, 20 x 6 in. (51 x 15cm)

Urological Drape



Sterile No. D1071 **Nonsterile No.** D1071NS **Qty. Cs**

- Adhesive: 24 x 2.7 in. (61 x 7cm)
- Aperture: 1 in. dia. (2.5cm dia.)
- Finger Cot Adhesive: 4 in. dia. (10cm dia.)

Cystoscopy Packs

Nonsterile No.	Description	29480 Cystoscopy Pack	29360 Cystoscopy Pack
29437NA	Cysto in.T in. drape, absorbent reinforcement	1	
29451N	Urological drape with absorbent reinforcement, attached leggings		1
8694-60	Suture bag	1	1
8520N	Absorbent towel	1	1
8421N	Pair of leggings, X-large	1	
8145N	Outer wrap/reinforced poly table cover 44 x 75 in. (112 x 190cm), reinforcement 24 x 75 in. (61 x 190cm)	1	1
	Quantity per case	14	10



Tiburon® Surgical Drapes and Packs

General surgery

Tiburon® Surgical Drapes and Packs

General surgery

General surgery

Tiburon® Head, Side and Foot Drape

Sterile No. Nonsterile No.

N/A	9467N
N/A	9465N
N/A	9468N
99 in. (257cm)	
	9467N
	53 in. (134cm)
	38 in. (96.5cm)
9465N	
	77 in. (195.5cm)
	9468N
	69 in. (175cm)

- Included in pack 9106

Tiburon® Universal Pack

Sterile No. Nonsterile No.

N/A	29667N
N/A	29665N
N/A	29669N
112 in. (285cm)	
	29667N
	72 in. (183cm)
	44 in. (112cm)
29665N	
	75 in. (190.5cm)
	29669N
	77 in. (195.5cm)

- Absorbent reinforcement with adhesive
- Cord-holding straps
- Included in pack 29118

Tiburon® Universal Split Pack

Sterile No. Nonsterile No.

29352	29351N	14
29440	29440NA	10
120 in. (305cm)		
		77 in. (195.5cm)
		53 in. (135cm)
		99 in. (251cm)

- Absorbent reinforcement
- Cord-holding tabs
- Split, 6 x 46 in. (15 x 117cm) (29440NA)
- Instrument pad (29351N)
- Included in pack 29131

Tiburon® Laparotomy Drape

Sterile No. Nonsterile No. Qty. Cs

29410	29410N	12
106 in. (269cm)		
		122 in. (310cm)
		77 in. (195.5cm)

- Absorbent reinforcement, 30.5 x 24 in. (77.5 x 61cm)
- Armboard covers
- Velcro® Cord-Holding Tabs
- Fenestration, 4 x 12 in. (10 x 30.5cm)
- Instrument pad
- Included in packs 29136

Tiburon® Laparotomy Drape with Pouches

Sterile No. Nonsterile No. Qty. Cs

9430	9430N	12
106 in. (269cm)		
		123.5 in. (314cm)
		77 in. (195.5cm)

- Absorbent reinforcement, 14 x 38.5 in. (35.5 x 98cm)
- Armboard covers
- Dual clear closed end poly troughs with Velcro® Closures, 8 x 38.5 in (20 x 98cm)
- Fenestration, 4 x 12 in. (10 x 30.5cm)

Tiburon® Thyroid Drape

Sterile No. Nonsterile No. Qty. Cs

29422	29422N	10
77 in. (195.5cm)		
		122 in. (310cm)

- Absorbent reinforcement, 22 x 30 in. (56 x 76cm)
- Cord-holding tabs
- Fenestration, diamond, 5 x 5 in. (13 x 13cm)
- Instrument pad



Tiburon® Surgical Drapes and Packs


General surgery

General surgery

Tiburon® Chest/Breast Drape

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

29420 29420N 12




- Absorbent reinforcement, 22 x 30 in. (56 x 76cm)
- Armboard covers
- Velcro® Cord-Holding Tabs
- Fenestration, 14 x 11 in. (35.5 x 28cm)
- Instrument pad

Tiburon® Laparoscopic Cholecystectomy Drape with Pouches

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

9428 9428NA 8




- Absorbent reinforcement, 36 x 42 in. (91 x 107cm)
- Armboard covers
- Cord-holding tabs
- Fabric pouches with adhesive, 10 x 42 in. (25 x 107cm)
- Fenestration, offset 5 in., 11 x 14 in. (12.7cm, 28 x 35.5cm)
- Instrument pad

Tiburon® General Endoscopy Drape

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

9458 9458N 12




- Absorbent reinforcement, 16 x 38.5 in. (40.5 x 98cm)
- Armboard covers
- Dual clear open ended poly troughs with Velcro® Closure, 11.125 x 30 in. (28 x 76cm)
- Fenestration, 10 x 11 in. (25.5 x 28cm)
- Velcro® Cord-Holding Tabs

Tiburon® Bariatric Drape

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

94590 94590N 8




- Armboard covers
- Fabric instrument troughs, 28.5 x 24 in. (72 x 61cm)
- Fenestration, chevron, 18 x 18 in. (45.5 x 45.5cm)
- Instrument pad
- Velcro® Cord-Holding Tabs

Tiburon® Major Abdominal Drape with Pouches

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

9457 9457N 12



- Absorbent reinforcement, 16 x 38.5 in. (40.5 x 98cm)
- Armboard covers
- Dual clear pouches with Velcro® Closure, 11.125 x 30 in. (28 x 76cm)
- Fenestration, 10 x 11 in. (25.5 x 28cm)
- Velcro® Cord-Holding Tabs

General surgery

Universal Packs

Nonsterile No.	Description	9106 Universal Pack	29118 Universal Pack	29131 Universal Split Pack
29440NA	Tiburon® Split Drape with absorbent reinforcement			1
29351N	Tiburon® Top Drape, absorbent reinforcement,instrument pad			1
29669N	Tiburon® Foot Drape, 75 x 77 in. (190.5 x 195.5cm), adhesive, absorbent reinforcement		1	
29667N	Tiburon® Head Drape, 112 x 72 in. (285 x 183cm), adhesive, absorbent reinforcement		1	
29665N	Tiburon® Side Drape, 44 x 77 in. (112 x 195.5cm), adhesive, absorbent reinforcement		2	
8520N	Absorbent towel	1	1	
8554N	Drape towel, adhesive, nonabsorbent		4	4
8337N	Mayo stand cover, reinforced poly, 23 in. (58cm) wide	1	1	1
8694-60	Suture bag	1	1	1
8146N	Outer wrap, reinforced poly table cover, 44 x 90 in. (112 x 228.6cm)	1	1	
8143N	Outer wrap, reinforced poly table cover, 50 x 90 in.(127 x 228.6cm)			1
		Quantity per case	5	6
				4

Laparotomy Packs

Nonsterile No.	Description	29136 Laparotomy Pack
9358N	Tiburon® Half Drape, 40 x 57 in. (101 x 145cm)	1
29410N	Tiburon® Laparotomy Drape, absorbent reinforcement	1
8554N	Drape towel, adhesive, nonabsorbent	4
8520N	Absorbent towel	1
8337N	Mayo stand cover, reinforced poly, 23 in. (58cm) wide	1
8694-60	Suture bag	1
8143N	Outer wrap, reinforced poly table cover, 50 x 90 in. (127 x 228.6cm)	1
		Quantity per case
		5




Tiburon® Surgical Drapes and Packs

Neurosurgery

Tiburon® Surgical Drapes and Packs


Neurosurgery



Neurosurgery

Tiburon® Neurological Drape


Sterile No. 9448 **Nonsterile No.** 9448N **Qty. Cs** 8



- Absorbent reinforcement, 18.5 x 28.5 in. (47 x 72cm)
- Cord-holding tabs
- Fenestration, oval with full incise, 11 x 7 in. (28 x 18cm)
- Fluid control drawstring pouch with screen and port, 20.5 x 40 in. (52 x 101.5cm)

Tiburon® Universal Spine Drape


Sterile No. 29418 **Nonsterile No.** 29418N **Qty. Cs** 7



- Absorbent reinforcement, 20 x 37 in. (51 x 94cm)
- Armboard covers
- Fenestration, 4 x 12 in. (10 x 30.5cm)
- Instrument pouches, 7.5 x 12 in. (19 x 30.5cm)
- Velcro® Cord-Holding Tabs

Tiburon® Craniotomy Split Drape

Sterile No. 9439 **Nonsterile No.** 9439N **Qty. Cs** 8



- Absorbent reinforcement, 20 x 38 in. (51 x 96.5cm)
- Attached overhead table cover, 80 x 90 in. (203 x 228.5cm)
- Cord-holding tabs
- Fenestration, diamond, 8 x 8 in. (20 x 20cm)
- Split, 54 in. (137cm)




Tiburon® Surgical Drapes and Packs

Obstetrics

Tiburon® Surgical Drapes and Packs

Obstetrics



Obstetrics

Tiburon® Cesarean Birth Drape

Sterile No.	Nonsterile No.	Qty. Cs
29449	29449N	7

106 in. (269cm)
121 in. (307cm)
77 in. (195.5cm)

- 360-degree fluid control pouch with contoured foam ring, Velcro® Cord-Holding Tab and ports
- Armboard covers
- Fenestration, circular with 2 in. (5cm) incise aperture, 16 in. (40.5cm)
- Reinforcement, 17 x 22 in. (43 x 56cm)
- Cord-holding tabs

Tiburon® Cesarean Birth Drape

Sterile No.	Nonsterile No.	Qty. Cs
29443	29443N	7

106 in. (269cm)
123.5 in. (313.7cm)
77 in. (195.5cm)

- 360-degree fluid control pouch with contoured foam ring, Velcro® Cord-Holding Tab and ports
- Armboard covers
- Fenestration, full incise, 9 x 14 in. (23 x 35.5cm)
- Reinforcement, 17 x 22 in. (43 x 56cm)
- Velcro® Cord-Holding Tabs

Tiburon® Cesarean Birth Drape

Sterile No.	Nonsterile No.	Qty. Cs
29448	29448N	7

106 in. (269cm)
123.5 in. (313.7cm)
77 in. (195.5cm)

- 360-degree fluid control pouch with contoured foam ring, Velcro® Cord-Holding Tab and ports
- Armboard covers
- Cord-holding tabs
- Fenestration, pear-shaped with full incise, 12 x 14 in. (30.5 x 35.5cm)
- Reinforcement, 17 x 22 in. (43 x 56cm)

Tiburon® Cesarean Birth Drape

Sterile No.	Nonsterile No.	Qty. Cs
29452	29452N	8

106 in. (269cm)
121 in. (307cm)
77 in. (195.5cm)

- 270-degree fluid control pouch with malleable wire
- Armboard covers
- Velcro® Cord-Holding Tabs
- Fenestration, pear-shaped with 2 in. (5cm) incise aperture, 12 x 14 in. (30.5 x 35.5cm)
- Reinforcement around fenestration (20.375 x 27 in. (51.8 x 68.6cm)) and below pouch (17 x 22 in. (43 x 56cm))

Tiburon® Cesarean Birth Drape

Sterile No.	Nonsterile No.	Qty. Cs
29447	29447N	8

106 in. (269cm)
128 in. (325cm)
77 in. (195.5cm)

- 270-degree fluid control pouch with malleable wire and ports
- Armboard covers
- Cord-holding tabs
- Fenestration, pear-shaped with full incise, 12 x 14 in. (30.5 x 35.5cm)
- Reinforcement around fenestration (20.375 x 27 in. (51.8 x 68.6cm)) and below pouch (17 x 22 in. (43 x 56cm))

Tiburon® Abdominal Cover

Sterile No.	Nonsterile No.	Qty. Cs
N/A	29357N	

57 in. (147cm)
39.5 in. (100cm)

- Absorbent reinforcement, 19 x 19 in. (48 x 48cm)
- Included in pack 29920

Obstetrics**Obstetrics Pack**


Nonsterile No.	Description	29920 Obstetrics Pack
29357N	Abdominal cover	1
8421N	Pair of leggings, X-large	1
8483N	Under-buttocks drape, reinforced poly	1
8365N	Setup cover, 44 x 71 in. (111.7 x 180cm)	1
8520N	Absorbent towel	2
8145N	Outer wrap/reinforced poly table cover 44 x 75 in. (111.7 x 190.5cm), reinforcement 24 x 75 in. (61 x 190.5cm)	1
	Quantity per case	6

Obstetrics



Tiburon® Surgical Drapes and Packs


Ophthalmology and ENT



Tiburon® Surgical Drapes and Packs
Ophthalmology and ENT

Ophthalmic Drape

Sterile No.	Nonsterile No.	Qty. Cs
7444	7444N	10




53 in. (134.5cm)

68 in. (178cm)

Ophthalmic Drape

Sterile No.	Nonsterile No.	Qty. Cs
7445	7445N	10




53 in. (134.5cm)

68 in. (178cm)

Ophthalmic Drape with Pouch

Sterile No.	Nonsterile No.	Qty. Cs
9562	9562N	10



77 in. (195.5cm)

110 in. (279cm)

Ophthalmology and ENT


- Cord-holding tab
- Dual saddlebag fluid control pouches with malleable strips
- Full incise, 4 x 5 in. (10 x 13cm)
- Malleable breather-bar
- SMS material

- Dual saddlebag fluid control pouches with malleable strips
- Malleable breather-bar
- Oval incise aperture, 1.5 x 2.5 in. (4 x 6cm)
- Cord-holding tab
- SMS material

- Attached fluid control saddlebag pouch
- Cord-holding tab
- Malleable breathing bar
- Oval aperture with incise, 1.5 x 2.5 in. (4 x 6cm)
- Separate incise sheet, 5 x 4 in. (13 x 10cm)
- SMS material

Ophthalmic Incise Drape

Sterile No.	Nonsterile No.	Qty. Cs
D2258	N/A	10




51 in. (129cm)

75 in.
(190cm)

Ophthalmic Aperture Drape

Sterile No.	Nonsterile No.	Qty. Cs
D1031	D1031NS	10




48 in. (121cm)

51 in.
(129cm)

Small Drape with Adhesive Aperture

Sterile No.	Nonsterile No.	Qty. Cs
8441	8441N	40



47 in. (119cm)


56 in. (142cm)

- Adhesive Area: 4.5 x 5 in. (11 x 12cm)
- Pouch: 21 x 11.5 in. (53 x 29cm)

- Aperture: 3 x 2 in. (7 x 5cm)
- Pouch: 11 x 12 in. (28 x 30cm)

- Oval fenestration
2 x 3 in. (5 x 7.5cm)
- SMS material

Ophthalmology and ENT




Ophthalmic Incise Drape

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

D1021	D1021NS	10
-------	---------	----

16 in. (40cm)




Ophthalmic Incise Drape

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

D1061	D1061NS	10
-------	---------	----

48 in. (121cm)




Ophthalmic Incise Pouch

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

D1025	D1025NS	10
-------	---------	----

19 in. (48cm)



- Adhesive Area: 5 x 4.25 in. (12 x 11cm)

- Adhesive Area: 5 x 4.25 in. (12 x 11cm)
- Pouch: 11 x 12 in. (28 x 30cm)


- Adhesive Area: 5 x 4.25 in. (12 x 11cm)
- Adhesive Tabs: (2) 1 x 2 in. (2.5 x 5cm)

Ophthalmic Incise Drape

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

D1060	D1060NS	10
-------	---------	----

48 in. (121cm)




Ophthalmic Incise Drape

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

D2019	D2019NS	30
-------	---------	----

42 in. (106cm)




Aperture Drape

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

D1020	D1020NS	10
-------	---------	----


16 in. (40cm)



- Adhesive Area: 5 x 4.25 in. (12 x 11cm)

- Adhesive Area: 3 in. dia. (7cm dia.)

- Aperture: 2.5 in. dia. (6cm dia.)




Tiburon® Surgical Drapes and Packs Ophthalmology and ENT

Ophthalmology and ENT


Refractive Drape

Sterile No.	Nonsterile No.	Qty. Cs
D1122		10




Ophthalmic Aperture Drape

Sterile No.	Nonsterile No.	Qty. Cs
D1024	D1024NS	10



Strabismus Drape

Sterile No.	Nonsterile No.	Qty. Cs
D2043	D2043NS	10




- Aperture: 1.5 x 1 in.
(3.8 x 2.5cm)

- Aperture: 2 x 1.5 in. (5cm x 4cm)
- Adhesive Tabs: (2) 1 x 2 in.
(2.5 x 5cm)

- Aperture: 5.2 x 1.5 in. (13 x 4cm)


Ophthalmic Aperture Drape

Sterile No.	Nonsterile No.	Qty. Cs
D2020	D2020NS	30



Tiburon® ENT Split Drape

Sterile No.	Nonsterile No.	Qty. Cs
N/A	9443N	




- Aperture: 3 x 2 in. (7 x 5cm)

- Split with adhesive, 2.5 x 25 in.
(6 x 63.5cm)
- Cord-holding tabs
- Instrument pad
- Included in pack 9195

Tiburon® Large ENT Split Drape

Sterile No.	Nonsterile No.	Qty. Cs
9447	9447N	8



- Instrument pad
- Split with adhesive, 2.5 x 25 in.
(6 x 63.5cm)
- Velcro® Cord-Holding Tabs
- Included in pack 9196

Tiburon® Surgical Drapes and Packs


Ophthalmology and ENT



Ophthalmology and ENT


Tiburon® Large ENT Split Drape

Sterile No.	Nonsterile No.	Qty. Cs
TB9441	9441N	8



Tiburon® Head/Bar Drape


Sterile No.	Nonsterile No.	Qty. Cs
9556	9556N	20



- Cord-holding tabs
- Turban
- Included in pack 9194

Tiburon® Bar Drape

Sterile No.	Nonsterile No.	Qty. Cs
N/A	9557N	



- Adhesive
- Included in pack 9195

- Absorbent reinforcement, 15.25 x 28 in. (39 x 71cm)
- Cord-holding tabs
- Instrument pad, 15.25 x 9.5 in. (38 x 24cm)
- Split with adhesive, 3.5 x 28 in. (9 x 71cm)

Ophthalmology and ENT

Head and Neck Packs

Nonsterile No.	Description	9195 Head and Neck Pack	9196 Head and Neck Pack
9556N	Tiburon® Head/Bar Drape, adhesive	1	
9447N	Tiburon® Large ENT split drape		1
9557N	Tiburon® Bar Drape, adhesive	1	
9443N	Tiburon® Split Drape, unreinforced, adhesive	1	
90270NA	Astound® Standard Scrub Nurse Gown with towel in outer wrap, large		1
9515NA	Astound® Standard Gown, large	1	1
8554N	Drape towel, adhesive, nonabsorbent	2	2
8520N	Absorbent towel		1
8337N	Mayo stand cover, reinforced poly, 23 in. wide (58cm)	1	1
8694-60	Suture bag	1	1
8146N	Outer wrap, reinforced poly table cover, 44 x 90 in. (111.7 x 228.6cm)	1	1
Quantity per case		8	6



Tiburon® Surgical Drapes and Packs

Orthopedics




Tiburon® Surgical Drapes and Packs

Orthopedics - Heavy Duty

Tiburon® HD Extremity Drape

Sterile No. 59405 **Nonsterile No.** 59405N **Qty. Cs** 6

117 in. (297cm)




- Absorbent reinforcement, 30.5 x 55.5 in. (77.5 x 141cm)
- Anesthesia screen
- Armboard covers
- Cord-holding tabs
- Fenestration, circular elastic, 2.5 in. (6cm)
- Heavyweight, premium fabric
- Included in pack 59180

Tiburon® HD Arthroscopy Drape

Sterile No. 59413 **Nonsterile No.** 59413N **Qty. Cs** 6

113 in. (287cm)




- Absorbent reinforcement
- Anesthesia screen
- Armboard covers
- Cord-holding tabs
- Dual fenestration, circular elastic, 2 in. (6cm)
- Fluid control pouch with malleable wire drawstring pouch and port
- Heavyweight, premium fabric
- Included in pack 9197

Tiburon® HD Hip Drape with Pouches

Sterile No. 59439 **Nonsterile No.** 59439N **Qty. Cs** 4

117 in. (297cm)



- Absorbent reinforcement, 21 x 57 in. (53 x 145cm)
- Anesthesia screen
- Armboard covers
- Heavyweight, premium fabric
- Clear leg positioning pouches, 54 x 30 in. (137 x 76cm)
- Fenestration, pear-shaped elastic, 6.75 x 8.25 in. (17 x 21cm)
- Velcro® Cord-Holding Tabs

Orthopedics - Heavy Duty

Tiburon® Surgical Drapes and Packs


Orthopedics



Orthopedics

Tiburon® Split Drape


Sterile No.	Nonsterile No.	Qty. Cs
29434	29434N	20



- Absorbent reinforcement, 20 x 18.5 in. (51 x 47cm)
- Split with adhesive, 1 x 28.25 in. (2.5 x 72cm)

Tiburon® Split Drape


Sterile No.	Nonsterile No.	Qty. Cs
29436	29436N	10



- Absorbent reinforcement, 34 x 30.5 in. (86 x 77.5cm)
- Split with adhesive, 2.5 x 37 in. (6 x 94cm)
- Included in packs 29182, 29183, 29184

Tiburon® Orthopedic Split Drape


Sterile No.	Nonsterile No.	Qty. Cs
29476	29476N	10



- Absorbent reinforcement, 30.5 x 42 in. (77.5 x 107cm)
- Split with adhesive, 6 x 54 in. (15 x 137cm)

Tiburon® Split Sheet Drape


Sterile No.	Nonsterile No.	Qty. Cs
29445	29445N	8



- Absorbent reinforcement, 33.5 x 53 in. (85 x 134.5cm)
- Split with adhesive, 4 x 40 in. (10 x 101cm)

Tiburon® Split Sheet Drape


Sterile No.	Nonsterile No.	Qty. Cs
29446	29446N	8



- Absorbent reinforcement, 33.5 x 53 in. (85 x 134.5cm)
- Split with adhesive, 6 x 46 in. (15 x 117cm)

Tiburon® Arthroscopy Drape for Dropped Leg Procedures

Sterile No.	Nonsterile No.	Qty. Cs
9414	9414NA	12




- Absorbent reinforcement, 20 x 20 in. (51 x 51cm)
- Dual fenestration, circular elastic, 2.5 in. (6cm)
- Fluid control pouch with malleable wire and port
- Cord-holding tabs



Tiburon® Surgical Drapes and Packs Orthopedics

Tiburon® Hip Drape


Sterile No. 29430 **Nonsterile No.** 29430N **Qty. Cs** 12



- Absorbent reinforcement, 30.5 x 50 in. (77.5 x 127cm)
- Cord-holding tabs
- Fenestration, circular elastic, 8 in. (20cm)
- Instrument pad
- Included in pack 29187

Tiburon® Hip Drape with Pouches


Sterile No. 29439 **Nonsterile No.** 29439N **Qty. Cs** 5



- Absorbent reinforcement, 21.375 x 57 in. (54 x 145cm)
- Anesthesia screen
- Armboard covers
- Clear leg positioning pouches, 54.5 x 30 in. (138 x 76cm)
- Fenestration, pear-shaped elastic, 6.75 x 8.25 in. (17 x 21cm)
- Velcro® Cord-Holding Tabs

Tiburon® Hand Drape

Sterile No. 29427 **Nonsterile No.** 29427N **Qty. Cs** 18



- Absorbent reinforcement, 16.75 x 34 in. (42.5 x 89cm)
- Fenestration, circular elastic, 2 in. (5cm).
- Hand table extension
- Velcro® Cord-Holding Tabs

Orthopedics

General Ortho Packs

Nonsterile No.	Description	29181 General Orthopedic Pack	29182 Orthopedic Pack
29432N	Tiburon® Split Drape with absorbent reinforcement, adhesive	1	
29351N	Tiburon® Top Drape with absorbent reinforcement, adhesive	1	
9349N	Tiburon® Fan-Folded Drape Sheet, 52 x 76 in. (132 x 193cm)		1
29436N	Tiburon® split drape, adhesive		1
9040N	Astound® Gown, poly reinforced, X-large	2	
8520N	Absorbent towel	2	
8694-60	Suture bag	1	1
8337N	Mayo stand cover, reinforced poly, 23 in. (58cm) wide	1	1
8475N	U-Drape, poly, adhesive, 60 x 84 in. (152 x 213cm), 8 x 30 in. (20 x 76cm) split	1	
8476NA	U-Drape, poly, adhesive, 60 x 72 in. (152 x 183cm), 6 x 21 in. (15 x 53cm) split		1
1582N	Stockinette, 12 x 48 in. (30.5 x 122cm), X-large	1	
8554N	Drape towel, adhesive, nonabsorbent	4	4
8143N	Outer wrap, reinforced poly table cover, 50 x 90 in. (127 x 228.6cm)	1	1
Quantity per case		4	5

Tiburon® Surgical Drapes and Packs


Orthopedics



Orthopedics

Tiburon® Extremity Drape


Sterile No. 29405 **Nonsterile No.** 29405N **Qty. Cs** 6



- Absorbent reinforcement, 30.5 x 45 in. (77.5 x 114cm)
- Anesthesia screen
- Armboard covers
- Cord-holding tabs
- Fenestration, circular elastic, 2.5 in. (6cm)

Tiburon® Extremity Drape


Sterile No. 29414 **Nonsterile No.** 29414N **Qty. Cs** 12



- Absorbent reinforcement, 30.5 x 45 in. (77.5 x 114cm)
- Cord-holding tabs
- Fenestration, circular elastic, 2.5 in. (6in.)
- Included in packs 29188, 29189

Tiburon® Bilateral Limb Drape


Sterile No. 29416 **Nonsterile No.** 29416N **Qty. Cs** 12



- Absorbent, reinforcement, 31 x 27 in. (79 x 68.5cm)
- Cord-holding tabs
- Dual fenestration, circular elastic, 1.5 in. (4cm)

Tiburon® Bilateral Limb Drape

Sterile No. 29417 **Nonsterile No.** 29417N **Qty. Cs** 8



- Absorbent, reinforcement, 37.5 x 61 in. (95 x 155cm)
- Anesthesia screen
- Armboard covers
- Cord-holding tabs
- Dual fenestration, circular elastic, 2.5 in. (6cm)

Tiburon® Beach Chair Shoulder Drape


Sterile No. 29369 **Nonsterile No.** 29369NA **Qty. Cs** 5




- Fenestration, pear-shaped, 5.5 x 5.75 in. (14 x 14.5cm)
- Fluid control pouch with malleable wire drawstring pouch and port
- Velcro® Cord-Holding Tabs

Tiburon® Shoulder Drape

Sterile No. N/A **Nonsterile No.** 9360N **Qty. Cs**



- Fenestration, pear-shaped with split and incise aperture, 11.5 x 14 in. (29 x 35.5cm)
- Fluid control pouch with drawstring and ports
- Cord-holding tabs
- Included in packs 29367




Tiburon® Surgical Drapes and Packs Orthopedics

Orthopedics

Tiburon® Split Drape

Sterile No.	Nonsterile No.	Qty. Cs
N/A	9361N	

77 in. (195.5cm)




116 in. (294.5cm)

- Reinforcement, 44 x 40 in. (112 x 101.5cm)
- Split with adhesive, 10 x 25 in. (25 x 63.5cm)
- Included in packs 9183, 29367, 29365, 29160

Tiburon® Top Drape

Sterile No.	Nonsterile No.	Qty. Cs
29352	29351N	14

99 in. (251.5cm)



53 in. (134.5cm)

- Absorbent reinforcement, 20 x 9.375 in. (51 x 24cm)
- Cord-holding tabs
- Instrument pad
- Included in packs 29181, 29131

Extremity Packs

Nonsterile No.	Description	29188 Fluid Control Extremity Pack
29414N	Tiburon® Extremity Drape	1
8694-60	Suture bag	1
8337N	Mayo stand cover, reinforced poly, 23 in. (58cm) wide	1
8375N	Plastic sheet, fan-folded, 44 x 60 in. (112 x 152cm)	1
1582N	Stockinette, 12 x 48 in. (30.5 x 122cm), X-large	
8554N	Drape towel, adhesive, nonabsorbent	2
8143N	Outer wrap, reinforced poly table cover, 50 x 90 in. (127 x 228.6cm)	1
Quantity per case		5

Shoulder Packs


Nonsterile No.	Description	29367 Fluid Control Shoulder Drape Pack	29365 Shoulder Drape Pack
9360N	Tiburon® Shoulder Drape with attached fluid control pouch, with incise	1	
9361N	Tiburon® Split Drape, poly reinforced, adhesive	1	1
		Quantity per case	10

Orthopedics



Patient Isolation Drape


Sterile No.	Nonsterile No.	Qty. Cs
D1014	D1014NS	5



- Adhesive Area: 9.5 x 19.5 in. (24 x 50cm)
- Adhesive Strip: 125 x 0.5 in. (317 x 1.3cm)

Patient Isolation Drape

Sterile No.	Nonsterile No.	Qty. Cs
D1017	D1017NS	5



- Adhesive Area 9.5 x 19.5 in. (24 x 50cm)
- Adhesive Strip 125 x 0.5 in. (317 x 1.3cm)
- Fluid Pouch 20 x 19 in. (51 x 48cm)
- Instrument Pouches (2) 6.7 x 11.5 in. (17 x 29cm)

PRODUCT KEY: Antimicrobial Adhesive



Adhesive




Matte Film




Clear Film




Tiburon® Surgical Drapes and Packs**Pacemaker****Fluoroscope Procedure Sheet**

Sterile No.	Nonsterile No.	Qty. Cs
D2221	D2221NS	20

**Pacemaker Drape**

Sterile No.	Nonsterile No.	Qty. Cs
29013	29013N	8



- Adhesive Area:
19.75 x 14 in. (43 x 33cm)
- Pouches:
11 x 23 in. (27 x 58cm)
8.5 x 11 in. (21 x 25cm)

- Absorbent reinforcement, 33.5 x 115 in.
(85 x 292cm)
- Clear side panels, 23 x 115 in.
(58 x 292cm)
- Fenestration, dual square, 6 in. (15cm),
with square incise aperture, 4 in. (10cm)



Tiburon® Surgical Drapes and Packs

Pediatrics

Tiburon® Surgical Drapes and Packs


Pediatrics



Pediatrics


Tiburon® Pediatric Drape

Sterile No.	Nonsterile No.	Qty. Cs
29492	29492N	10



Tiburon® Neonatal Drape


Sterile No.	Nonsterile No.	Qty. Cs
29494	29494N	20



Tiburon® Split Drape

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

29434	29434N	20
-------	--------	----



- Absorbent reinforcement, 18 x 20 in. (46 x 51cm)
- Fenestration, 3 x 4 in. (7.5 x 10cm)

- Absorbent reinforcement, 20.375 x 18.5 in. (52 x 47cm)
- Split with adhesive, 1 x 28.25 in. (2.5 x 72cm)

- Absorbent reinforcement, 19 x 20.375 in. (48 x 52cm)
- Cord-holding tabs
- Fenestration, 6 x 6 in. (15 x 15cm)
- Instrument pad

Tiburon® Pediatric Directional Incision Drape

Sterile No.	Nonsterile No.	Qty. Cs
-------------	----------------	---------

9777	9777N	60
------	-------	----



- Fenestration, circular, 4 in. (10cm)



Tiburon® Surgical Drapes and Packs

Peri-GYN

Tiburon® Surgical Drapes and Packs


Peri-GYN



Peri-GYN

Tiburon® Lithotomy Drape


Sterile No.	Nonsterile No.	Qty. Cs
29455	29455N	8



- Absorbent reinforcement, 20.375 x 19 in.
- Anesthesia screen
- Attached leggings
- Cord-holding tabs
- Fluid control pouch with filtering screen and suction port
- Rectangular perineal fenestration with adhesive, 2.5 x 6.5 in. (6 x 16.5cm)

Tiburon® LAVH Drape

Sterile No.	Nonsterile No.	Qty. Cs
29474	29474N	8



- Absorbent reinforcement, 12 x 19.25 x 16.75 in. (30.5 x 49 x 42.5cm)
- Anesthesia screen
- Armboard covers
- Attached leggings with tieback straps
- Fluid control pouch with closed suction port and drawstring
- Instrument pad
- Perineal fenestration with adhesive, 3 x 6 in. (7.5 x 15cm)
- Perineal flap
- Trapezoidal abdominal fenestration with adhesive, 7 x 9 x 9 in. (18 x 23 x 23cm)
- Velcro® Cord-Holding Tabs
- Included in pack 29475

Peri-GYN Packs


Nonsterile No.	Description	29242 Laparoscopy Pack	29146 Abdominal Lithotomy Pack	29276 Peri-GYN Pack	9165 Peri-GYN Pack
29438N	Tiburon® Lithotomy "T" Drape with absorbent reinforcement, adhesive			1	
29402N	Tiburon® Laparoscopy Drape with absorbent reinforcement, adhesive	1			
29300N	Abdominal lithotomy drape		1		
9359N	Drape sheet, adhesive				1
8483N	Under-buttocks drape, reinforced poly		1		1
8554N	Drape towel, adhesive, nonabsorbent	1		1	1
8520N	Absorbent towel	1		1	1
8421N	Pair of leggings, X-large	1		1	1
8694-60	Suture bag	1		1	1
8146N	Outer wrap, reinforced poly table cover, 44 x 90 in. (111.8 x 228cm)	1		1	1
	Quantity per case	10	8	10	12



Tiburon® Surgical Drapes and Packs


Special procedures

Special procedures



Special Procedure Sheet


Sterile No.	Nonsterile No.	Qty. Cs
D2224	D2224NS	20



- Adhesive Areas: 8.25 x 10 in. (20 x 15cm)
- Aperture: 3.5 x 4.5 in (8 cm x 11 cm)

Fistula Drape

Sterile No.	Nonsterile No.	Qty. Cs
D2264	D2264NS	10



- Aperture: 3 x 9 in. (8 x 23cm)



Tiburon® Surgical Drapes and Packs

Standard packs

Standard packs

Universal Packs

Universal Pack

Sterile No. 9106

NS Cat. No.	Contents	Qty. Cs
9465N	2 Tiburon® Side Drape, adhesive	5
9468N	1 Tiburon® Foot Drape, adhesive	
9467N	1 Tiburon® Head Drape, adhesive	
8520N	1 Absorbent towel	
8337N	1 Mayo stand cover, reinforced poly, 23 in. (58cm) wide	
8694-60	1 Suture bag	
8146N	1 Outer wrap, reinforced poly table cover, 44 x 90 in. (111.8 x 228.6cm)	

Universal Pack

Sterile No. 29118

NS Cat. No.	Contents	Qty. Cs
29669N	1 Tiburon® Foot Drape with adhesive, 77 x 75 in. (195.6 x 190.5cm)	6
29667N	1 Tiburon® Head Drape with adhesive, 72 x 112 in. (182.9 x 284.5cm)	
29665N	2 Tiburon® Side Drape with tape, 44 x 77 in. (111.8 x 195.6cm)	
8554N	4 Drape towel, adhesive, nonabsorbent	
8337N	1 Mayo stand cover, reinforced poly, 23 in. (58cm) wide	
8694-60	1 Suture bag	
8146N	1 Outer wrap, reinforced poly table cover, 44 x 90 in. (111.8 x 228.6cm)	
8520N	1 Absorbent towel	

Universal Split Pack

Sterile No. 29131

NS Cat. No.	Contents	Qty. Cs
29440NA	1 Tiburon® Split Drape with absorbent reinforcement	4
29351N	1 Tiburon® Top Drape, absorbent reinforcement/instrument pad 53 x 99 in. (134.6 x 251.5cm)	
8554N	4 Drape towel, adhesive, nonabsorbent	
8337N	1 Mayo stand cover, reinforced poly, 23 in. (58cm) wide	
8694-60	1 Suture bag	
8143N	1 Outer wrap, reinforced poly table cover, 50 x 90 in. (127 x 228.6cm)	

Laparotomy Packs

Laparotomy Pack

Sterile No. 29136

NS Cat. No.	Contents	Qty. Cs	5
9358N	1 Tiburon® Fan-Folded Drape Sheet, 40 x 57 in. (101 x 145cm)		
29410N	1 Tiburon® Laparotomy Drape with absorbent reinforcement and armboard covers		
8554N	4 Drape towel, adhesive, nonabsorbent		
8520N	1 Absorbent towel		
8337N	1 Mayo stand cover, reinforced poly, 23 in. (58cm) wide		
8694-60	1 Suture bag		
8143N	1 Outer wrap, reinforced poly table cover, 50 x 90 in. (127 x 228.6cm)		

Cystoscopy Packs

Cystoscopy Pack

Sterile No. 29480

NS Cat. No.	Contents	Qty. Cs	14
8694-60	1 Suture bag		
8520N	1 Absorbent towel		
8421N	1 Pair of leggings, X-large		
29437NA	1 Cystoscopy T-drape		
8145N	1 Outer wrap, reinforced poly table cover, 44 x 75 in. (111.8 x 190.5cm)		

Cystoscopy Pack

Sterile No. 29360

NS Cat. No.	Contents	Qty. Cs	10
8694-60	1 Suture bag		
8520N	1 Absorbent towel		
29451N	1 Tiburon® Urological Drape		
8145N	1 Outer wrap, reinforced poly table cover, 44 x 75 in. (111.8 x 190.5cm)		

Obstetrics Pack

Obstetrics Pack

Sterile No. 29920

NS Cat. No.	Contents	Qty. Cs	6
29357N	1 Tiburon® Abdominal Cover with absorbent reinforcement		
8421N	1 Pair of leggings, X-large		
8483N	1 Under-buttocks drape, reinforced poly		
8365N	1 Setup cover, 44 x 71 in. (111.8 x 180.3cm)		
8520N	2 Absorbent towel		
8145N	1 Outer wrap/reinforced poly table cover, 44 x 75 in. (111.8 x 190.5cm), reinforcement 24 x 75 in. (61 x 190.5cm)		

Standard packs

Head and Neck Packs

Head and Neck Pack

Sterile No. 9195

NS Cat. No.	Contents	Qty. Cs 8
9557N	1 Tiburon® Bar Drape, adhesive	
9443N	1 Tiburon® Split Sheet, unreinforced, adhesive	
9515NA	1 Astound® Standard Gown, sterile back	
8554N	2 Drape towel, adhesive, nonabsorbent	
8337N	1 Mayo stand cover, reinforced poly, 23 in. (58cm) wide	
8694-60	1 Suture bag	
8146N	1 Outer wrap, reinforced poly table cover, 44 x 90 in. (111.8 x 228.6cm)	

Head and Neck Pack

Sterile No. 9196

NS Cat. No.	Contents	Qty. Cs 6
9556N	1 Tiburon® Head/Bar Drape, adhesive	
9447N	1 Tiburon® Large ENT Split Drape	
90270NA	1 Astound® Standard Scrub Nurse Gown, large, with towel in outer wrap	
9515NA	1 Astound® Standard Gown, sterile back, large	
8520N	1 Absorbent towel	
8337N	1 Mayo stand cover, reinforced poly, 23 in. (58cm) wide	
8694-60	1 Suture bag	
8554N	2 Drape towel, adhesive, non-absorbant	
8146N	1 Outer wrap, reinforced poly table cover, 44 x 90 in. (111.8 x 228.6cm)	

Standard packs

Standard packs

Orthopedic Packs

General Orthopedic Pack

Sterile No. 29181

NS Cat. No.	Contents	Qty. Cs 4
29432N	1 Tiburon® Split Drape with absorbent reinforcement, adhesive	
29351N	1 Tiburon® Top Drape with absorbent reinforcement 53 x 99 in. (134.6 x 251.5cm), adhesive	
9040N	2 Astound® Reinforced Gown, X-large	
8520N	2 Absorbent towel	
8694-60	1 Suture bag	
8337N	1 Mayo stand cover, reinforced poly, 23 in. (58cm) wide	
8475N	1 U-Drape, poly, adhesive, 60 x 84 in. (152.4 x 213.4cm)	
1582N	1 Stockinette, 12 x 48 in. (30.5 x 121.9cm), X-large	
8554N	4 Drape towel, adhesive, nonabsorbent	
8143N	1 Outer wrap, reinforced poly table cover, 50 x 90 in. (127 x 228.6cm)	

Orthopedic Pack

Sterile No. 29182

NS Cat. No.	Contents	Qty. Cs 5
9349N	1 Tiburon® Fan-Folded Drape Sheet, 52 x 76 in. (132 x 193cm)	
29436N	1 Tiburon® Split Drape, adhesive	
8694-60	1 Suture bag	
8337N	1 Mayo stand cover, reinforced poly, 23 in. (58cm) wide	
8476NA	1 Prep U-drape, 60 x 72 in. (152.4 x 182.9cm)	
A8554N	4 Drape towel, adhesive, nonabsorbent	
A8143N	1 Outer wrap, reinforced poly table cover, 50 x 90 in. (127 x 228.6cm)	

Standard packs

Extremity Packs

Extremity Pack

Sterile No. 29188

NS Cat. No.	Contents	Qty. Cs	5
29414N	1 Tiburon® Extremity Drape		
8694-60	1 Suture bag		
8337N	1 Mayo stand cover, reinforced poly, 23 in. (58cm) wide		
8375N	1 Plastic sheet, fan-folded, 44 x 60 in. (111.8 x 152.4cm)		
8554N	2 Drape towel, adhesive, nonabsorbent		
8143N	1 Outer wrap, reinforced poly table cover, 50 x 90 in. (127 x 228.6cm)		

Shoulder Packs

Fluid Control Shoulder Drape Pack

Sterile No. 29367

NS Cat. No.	Contents	Qty. Cs	10
9360N	1 Tiburon® Shoulder Drape with attached fluid control pouch, with incise		
9361N	1 Tiburon® Split Drape, poly reinforced, adhesive		

Laparoscopy Packs

Laparoscopy Pack

Sterile No. 29242

NS Cat. No.	Contents	Qty. Cs	10
29402N	1 Tiburon® Laparoscopy Drape with absorbent reinforcement, adhesive		
8554N	1 Drape towel, adhesive, nonabsorbent		
8520N	1 Absorbent towel		
8421N	1 Pair of leggings, X-large		
8694-60	1 Suture bag		
8146N	1 Outer wrap, reinforced poly table cover, 44 x 90 in. (111.8 x 228.6cm)		

Abdominal Lithotomy Pack

Sterile No. 29146

NS Cat. No.	Contents	Qty. Cs	8
29300N	1 Tiburon® Abdominal Lithotomy Drape		
8483N	1 Under-buttocks drape, reinforced poly		

Peri-GYN and Laparoscopy Packs**Peri-GYN Pack****Sterile No. 29276**

NS Cat. No.	Contents	Qty. Cs	10
29438N	1 Tiburon® Lithotomy "T" Drape with absorbent reinforcement, adhesive		
8554N	1 Drape towel, adhesive, nonabsorbent		
8520N	1 Absorbent towel		
8421N	1 Pair of leggings, X-large		
8694-60	1 Suture bag		
8146N	1 Outer wrap, reinforced poly table cover, 44 x 90 in. (111.8 x 228.6cm)		

Peri-GYN Pack**Sterile No. 9165**

NS Cat. No.	Contents	Qty. Cs	12
9359N	1 Drape sheet, adhesive 38.5 x 56 in. (97.8 x 142.2cm)		
8483N	1 Under-buttocks drape, reinforced poly		
8421N	1 Pair of leggings, X-large		
8520N	1 Absorbent towel		
8554N	1 Non-absorbent towel with adhesive strip		
8694-60	1 Suture bag		
8146N	1 Outer wrap, reinforced poly table cover, 44 x 90 in. (111.8 x 228.6cm)		

Sterile Number Index

1584.....	10	29476.....	40
1585.....	10	29480	22, 55
1586.....	10	29492	48
1587.....	10	29494	48
29013.....	46	29527	20
29118.....	26, 54	29529	14
29131	26, 54	29570	14
29136.....	26, 55	29920	31, 55
29146.....	50, 58	32474	5
29153.....	18	3524.....	11
29154.....	17	39015	5
29181.....	41, 57	39045	5
29182.....	41, 57	39049	5
29188.....	43, 58	39075	5
29242.....	50, 58	39079	5
29276.....	50, 59	39099	5
29322.....	20	59405	39
29352.....	24, 43	59413	39
29360.....	22, 55	59439	39
29365.....	43	599	5
29367.....	43, 58	7444.....	33
29369.....	42	7445.....	33
29405.....	42	8186.....	12
29410.....	24	8337.....	12
29414.....	42	8339.....	12
29416.....	42	8376.....	12
29417.....	42	8377.....	12
29418.....	28	8378.....	12
29420.....	25	8420.....	10
29422.....	24	8430.....	10
29427.....	41	8441.....	33
29430.....	41	8475.....	10
29434.....	40, 48	8476.....	10
29436.....	40	8482.....	9
29439.....	41	8483.....	9
29440.....	24	8484.....	9
29443.....	30	8487.....	9
29445.....	40	89005	5
29446.....	40	89015	5
29447.....	30	89045	5
29448.....	30	89075	5
29449.....	30	9010.....	5
29452.....	30	9040.....	5
29455.....	22, 50	9070.....	5
29456.....	14	9106.....	26, 54
29458.....	15	9154.....	17
29460.....	14	9155.....	17
29474.....	50	9156.....	17

Sterile Number Index

Sterile Number Index

9162.....	17	D1010	8
9165.....	50, 59	D1012	12
9195.....	37, 56	D1014	44
9196.....	37, 56	D1015	10
9349.....	8	D1016	11
9355.....	8	D1017	44
9358.....	8	D1018	11
9414.....	40	D1018L.....	11
9428.....	25	D1019	10
9430.....	24	D1020	34
9439.....	28	D1021	34
9444.....	8	D1024	35
9446.....	8	D1025	34
9447.....	35	D1031	33
9448.....	28	D1060	34
9457.....	25	D1061	34
9458.....	25	D1071	22
94590.....	25	D1072	11
9461.....	9	D1089	8
9462.....	8	D1122	35
9556.....	36	D1195	14
9562.....	33	D2008A	12
9777.....	48	D2019	34
		D2020	35
		D2043	35
		D2158	12
		D2221	46
		D2224	52
		D2244	14
		D2258	33
		D2264	52

A

ASG9511	5
ASG9515	5
ASG9541	5
ASG9541EL.....	5
ASG9545	5
ASG9575	5

D

D1000	8
D1003	11

T

TB9441.....	36
-------------	----

Nonsterile Number Index

1582N	10, 41, 43, 57	29492N.....	48
1585N	10	29494N.....	48
1586N	10	29527N.....	20
29013N.....	46	29529N.....	14
29154N.....	17	29570N.....	14
29300N.....	50, 58	29665N.....	24, 26, 54
29322N.....	20	29667N.....	24, 26, 54
29351N.....	24, 26, 41, 43, 54, 57	29669N.....	24, 26, 54
29357N.....	30, 31, 55	32474N.....	5
29365.....	43	39015N.....	5
29369NA.....	42	39045N.....	5
29402N.....	50, 58	39049N.....	5
29405N.....	42	39075N.....	5
29410N.....	24, 26, 55	39079N.....	5
29414N.....	42-43, 58	39099N.....	5
29416N.....	42	59405N.....	39
29417N.....	42	59413N.....	39
29418N.....	28	59439N.....	39
29420N.....	25	7444N.....	33
29422N.....	24	7445N.....	33
29427N.....	41	7554N.....	18
29430N.....	41	8143N.....	18, 26, 41, 43, 54-55, 57-58
29432N.....	41, 57	8145N.....	12, 22, 31, 55
29434N.....	40, 48	8146N.....	12, 26, 37, 50, 54, 56, 58-59
29436N.....	40-41, 57	8152N.....	18
29437NA.....	22, 55	8184N.....	12
29438N.....	50, 59	8186N.....	12
29439N.....	41	8188N.....	18
29440NA.....	24, 26, 54	8337N.....	12, 18, 26, 37, 41, 43, 54-58
29443N.....	30	8339NI.....	12
29445N.....	40	8365N.....	31, 55
29446N.....	40	8375N.....	43, 58
29447N.....	30	8420NA.....	10
29448N.....	30	8421N.....	22, 31, 50, 55, 58-59
29449N.....	30	8430NA.....	10
29451N.....	22, 55	8441N.....	33
29452N.....	30	8470N.....	18
29455N.....	50	8475N.....	10, 41, 57
29455NA.....	22	8476NA.....	10, 41, 57
29456N.....	14	8482N.....	9
29458N.....	15	8483N.....	9, 31, 50, 55, 58, 59
29460N.....	14	8484N.....	9
29462N.....	18	8487N.....	9
29464N.....	18	8520N.....	18, 22, 26, 31, 37, 41, 50, 54-59
29467N.....	18	8554N.....	18, 26, 37, 41, 43, 50, 54-59
29468N.....	18	8694-60.....	18, 22, 26, 37, 41, 43, 50, 54-59
29474N.....	50	89005N.....	5
29476N.....	40	89015N.....	5

Nonsterile Number Index

89045N	5
89075N	5
9010NA	5
90270NA	37, 56
9040N	18, 41, 57
9040NA	5
90680NA	18
9070N	5
90780N	18
9150N	18
9154N	17-18
9155N	17
9156N	17
9158N	18
9162N	17
9349N	8, 18, 41, 57
9355N	8
9358N	8, 26, 55
9359N	50, 59
9360N	42-43, 58
9361N	43, 58
9414NA	40
9428NA	25
9430N	24
9439N	28
9441N	36
9443N	35, 37, 56
9444N	8
9446N	8
9447N	35, 37, 56
9448N	28
9457N	25
9458N	25
94590N	25
9460N	18
9461N	9
9462N	8
9465N	24, 54
9467N	24, 54
9468N	24, 54
9511N	5
9515NA	5, 37, 56
9541ELNA	5
9541N	5, 18
9545NA	5
9556N	36-37, 56
9557N	36-37, 56
9562N	33
9575NA	5
9777N	48

A

A8143N	57
A8554N	57

D

D1000NS	8
D1003NS	11
D1010NS	8
D1012NS	12
D1014NS	44
D1015NS	10
D1016NS	11
D1017NS	44
D1018LNS	11
D1018NS	11
D1019NS	10
D1020NS	34
D1021NS	34
D1024NS	35
D1025NS	34
D1031NS	33
D1060NS	34
D1061NS	34
D1071NS	22
D1072NS	11
D1089NS	8
D1195NS	14
D2008NS	12
D2019NS	34
D2020NS	35
D2043NS	35
D2158NS	12
D2221NS	46
D2224NS	52
D2244NS	14
D2264NS	52

Product Index

B

Bags

Isolation Bags..... 11

C

Covers

Equipment Covers..... 12
 Heavy-Duty Back Table Covers 12
 Mayo Stand Covers 12
 Standard Back Tables Covers..... 12

D

Drapes

Angiography

Angiography Sheet..... 14
 Extra Large Femoral Angiography Drape 14
 Radial/Femoral Angiography Drape..... 14
 Tiburon® Brachial Angiography Drape 15
 Tiburon® Femoral Angiography Drape..... 14

Cardiovascular

Tiburon® Bottom Drape with Reinforcement..... 18
 Tiburon® Cardiovascular Incise Drape 17
 Tiburon® Cardiovascular Peri-Groin Drape..... 17
 Tiburon® Cardiovascular Split Drape..... 17
 Tiburon® Cardiovascular Split Drape with Pouches..... 17
 Tiburon® Head, Side and Foot Drapes 18

Central line

Full Body Central Line Drape..... 20
 Full Body Central Line Insertion Drape..... 20

Cystoscopy

Tiburon® Cystoscopy "T" Drape 22
 Tiburon® Lithotomy Drape with Fluid Control Pouch..... 22
 Tiburon® Urological Drape..... 22
 Urological Drape..... 22

General surgery

Tiburon® Bariatric Drape..... 25
 Tiburon® Chest/Breast Drape..... 25
 Tiburon® Craniotomy Split Drape..... 28
 Tiburon® General Endoscopy Drape..... 25
 Tiburon® Head, Side and Foot Drape 24
 Tiburon® Laparoscopic Cholecystectomy Drape with Pouches..... 25
 Tiburon® Laparotomy Drape 24
 Tiburon® Laparotomy Drape with Pouches 24
 Tiburon® Major Abdominal Drape with Pouches 25

Drapes, continued

General surgery

Tiburon® Neurological Drape..... 28
 Tiburon® Thyroid Drape 24
 Tiburon® Universal Pack..... 24
 Tiburon® Universal Spine Drape 28
 Tiburon® Universal Split Pack..... 24

Obstetrics

Tiburon® Abdominal Cover..... 30
 Tiburon® Cesarean Birth Drape 30

Ophthalmology and ENT

Aperture Drape 34
 Ophthalmic Aperture Drape 33, 35
 Ophthalmic Drape 33
 Ophthalmic Drape with Pouch..... 33
 Ophthalmic Incise Drape..... 33, 34
 Ophthalmic Incise Pouch 34
 Refractive Drape 35
 Small Drape with Adhesive Aperture..... 33
 Strabismus Drape 35
 Tiburon® Bar Drape 36
 Tiburon® ENT Split Drape 35
 Tiburon® Head/Bar Drape 36
 Tiburon® Large ENT Split Drape 35-36

Orthopedics

Patient Isolation Drape..... 44
 Tiburon® Arthroscopy Drape for Dropped Leg Procedures 40
 Tiburon® Beach Chair
 Shoulder Drape 42
 Tiburon® Bilateral Limb Drape 42
 Tiburon® Extremity Drape 42
 Tiburon® Hand Drape 41
 Tiburon® Hip Drape 41
 Tiburon® Hip Drape with Pouches 41
 Tiburon® Orthopedic Split Drape 40
 Tiburon® Shoulder Drape 42
 Tiburon® Split Drape 40, 43
 Tiburon® Split Sheet Drape 40
 Tiburon® Top Drape 43

Orthopedics - Heavy Duty

Tiburon® HD Arthroscopy Drape 39
 Tiburon® HD Extremity Drape 39
 Tiburon® HD Hip Drape with Pouches 39

Pacemaker

Fluoroscope Procedure Sheet..... 46
 Pacemaker Drape 46

Product Index

Drapes, continued

Pediatrics

Tiburon® Neonatal Drape	48
Tiburon® Pediatric Directional Incision Drape	48
Tiburon® Pediatric Drape	48
Tiburon® Split Drape	48

Peri-GYN

Tiburon® LAVH Drape.....	50
Tiburon® Lithotomy Drape.....	50

Poly U-Drapes	10
---------------------	----

Special procedures

Fistula Drape	52
Special Procedure Sheet.....	52

Under-Buttocks Drapes.....	9
----------------------------	---

G

Gowns, Surgical

Astound® Fabric-Reinforced Gowns	5
Astound® Poly-Reinforced Gowns.....	5
Astound® Standard Gowns	5

SmartGown™ Surgical Gowns.....	5
--------------------------------	---

L

Leggings

Leggings.....	10
---------------	----

P

Packs

Abdominal Lithotomy Packs	50, 58
Cardiovascular Packs.....	18
Cardiovascular Pack with reinforced Gowns	18
Cystoscopy Packs.....	22, 55

Packs, continued

Extremity Packs	43, 58
Fluid Control Extremity Packs.....	43
Fluid Control Shoulder Drape Packs.....	58
General Orthopedic Packs.....	41, 57
Head and Neck Packs.....	37, 56
Laparoscopy Packs.....	50, 58
Laparotomy Packs	26, 55
Lithotomy Packs.....	58
Obstetrics Packs	31, 55
Orthopedic Packs	41, 57
Peri-GYN Packs	50, 59
Shoulder Packs	43, 58
Universal Packs.....	26, 54
Universal Split Packs.....	26, 54

Pouches

Pouches.....	11
--------------	----

S

Sheets

Drape Sheets.....	8
Drape Sheet, Square-Folded	9

Sleeves

Sterile Disposable Sleeve.....	5
--------------------------------	---

Stockinettes

Stockinettes.....	10
-------------------	----

T

Tape

Adhesives and Drape Tape.....	11
-------------------------------	----

Towels

Drape Towels.....	8
-------------------	---

This literature piece is for distribution through ISKUS HEALTH in the United Kingdom and Ireland only.

© 2015 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE, ASTOUND, SMARTGOWN AND TIBURON are trademarks or registered trademarks of Cardinal Health.
All other marks are the property of their respective owners. Lit. No. 2PERI15-269739 (09/2015). PRC (09/2015).



Distributed by



Customer Services and Distribution Contact Details

Ireland

Tel: +353 (0) 1 428 7895/6
Fax: +353 (0) 1 428 7876
Email: info@iskushealth.com

United Kingdom

Tel: +44 (0) 208 869 6509
Fax: +44 (0) 208 869 6561
Email: uksales@iskushealth.com



www.iskushealth.com