



ISKUS HEALTH

IMPROVING PATIENT OUTCOMES

SILICONE WOUND DRAINS



For over 20 year Iskus Health has been supplying high-quality, reliable wound drainage products which perform to as our customers expect and hence we have become a trusted solution for these products.

Partnering with PRIMED GmbH who is a highly experienced manufacturer of high-quality wound drainage products gives us confidence that the products will perform perfectly for you and your patients with the reassurance of German medical engineering.

All products are developed and manufactured in-house employing a quality standard certified in accordance with DIN EN 13485. From with the initial product idea to material tests and tool construction right up to the product launch, PRIMED produce innovative medical devices from a single source. Production spans the entire manufacturing process from injection moulding/extrusion/assembly up to sterilisation.


Explanations

STERILE EO	Sterilised with ethylene oxide
PU/SU	Packaging unit / shipping unit
REF	Reference number/product Code
Ch..	French scale = Unit indicating the External diameter of drains (1 Ch. = 1/3 mm)
∅	Diameter

LOW-VACUUM SYSTEMS

Wound Drainage Kit with flat drain

Low-vacuum drains with flat drains are used for suction-assisted, post-operative drainage of wound secretions. They comprise a silicone collecting bulb, which can be compressed to generate a vacuum for active suction support and a J.P. flat drain made of silicone. In addition to this active drainage, the system can also function by means of passive (gravity) drainage.



Drainage kit with flat drain

- Silicone collecting bulb 100 cm³
 - Volume = 100 ml
 - Transparent material allows visual control of drained fluid
 - With universal hanger
 - With non-return valve
 - With scale to provide an at-a-glance information about current fill levels
 - Stepless plug connection for safe and gentle removal of the drain
- J.P. flat drain
 - Soft, medical grade silicone ensures good tissue compatibility
 - Total length = 600 mm, of which 200 mm are a fully perforated flat drain
 - Fully radiopaque
 - Inner ribbing prevents the drain from collapsing (even at higher pressures or when a vacuum is applied)
- With rounded, atraumatic holes


Flat drain			Drain body		PU/SU	REF
Width	Height	Length	Ch.	Length		
4 mm	3mm	200mm	Ch.15	600mm	10/40	24840
7 mm	3mm	200mm	Ch.15	600mm	10/40	24841
10 mm	4mm	200mm	Ch.15	600mm	10/40	24842



Silicone reservoir

- With connector for J.P. flat drains
- Made of medical grade silicone

PU/SU	REF	
100 cm ³	10/30	21078
150 cm ³	10/30	21085
400 cm ³	10/30	21086



Flat drain

LOW-VACUUM SYSTEMS

Silicone reservoir with Luer lock connection

This collection reservoir can be used in combination with a drain for suction-assisted, postoperative drainage of wound secretions. The Silicone collection reservoir functions by means of passive (gravity) drainage, but can also be compressed to generate a negative pressure and support active suction.



Product features


- Silicone collecting bulb 400 cm³
 - Volume = 400 ml
 - Transparent material allows visual control of drained fluid
 - With universal hanger
 - With non-return valve
 - With scale to provide an at-a-glance information about current fill levels
 - With Luer lock connection
 - Stepless plug connection for safe and gentle removal of silicone drains
- Venting filter

	PU/SU	REF
400 cm ³	15/15	21203

SILICONE ROUND DRAINS

Soft Silicone Round Drains


Soft drains are round drains made of medical grade silicone and used for open or half-open wound drainage via gravity.



Product features:

- Soft, medical grade silicone ensures good tissue compatibility
- Rounded, atraumatic tip
- 6 side openings
- With radiopaque stripe
- Length = 500 mm

Ch.	PU/SU	REF
8	25/100	21899
10	25/100	21900
12	25/100	21901
14	25/100	21130
15	25/100	21902
18	25/100	21903
20	25/100	21131
21	25/100	21904
24	25/100	21905
26	25/100	21132
27	25/100	21906
30	25/100	21907
33	25/100	21908
36	25/100	21909
39	25/100	21133




Soft drain

SILICONE CHANNEL DRAINS

without trocar

Channel drains without trocar


Cross-profile drain made of silicone with cross-shaped, slotted drain profile for open or half-open drainage of wound secretion.



Product features:

- Soft, medical grade silicone ensures good tissue compatibility
- With radiopaque stripe
- Kink resistant thanks to channel structure
- 4 drainage channels ensure continuous secretion flow and avoid incrustations
- Easy to remove thanks to low incidence of tissue ingrowth
- Length markings facilitate repositioning
- Transparent drain tube allows for easy visual inspection

Ch.	Total length	Length of drainage channels	PU/SU	REF
10	1100 mm	300 mm	10/10	24458
15	1100 mm	300 mm	10/10	24468
19	1100 mm	300 mm	10/10	24472
21	1100 mm	300 mm	10/10	24469
24	1100 mm	300 mm	10/10	24470
27	1100 mm	300 mm	10/10	24471




Channel drain

SILICONE CHANNEL DRAINS

with trocar

Channel drains with trocar


Channel drain made of silicone with cross-shaped, slotted drain profile and preconnected single-use guide needle (trocar) for optimal positioning. These drains are used for open and half-open drainage of wound secretion.



Product features:

- Made of soft medical grade silicone and a preconnected guiding needle (trocar) made of medical grade stainless steel
- Sharp trocar enables ideal positioning of the drain
- With radiopaque stripe
- Kink resistant thanks to channel structure
- 4 drainage channels ensure continuous secretion flow and avoid incrustations
- Easy to remove thanks to low incidence of tissue ingrowth
- Length markings facilitate repositioning
- Transparent drain tube allows for easy visual inspection

Ch.	Total length	Length of drainage channels	PU/SU	REF
10	1100 mm	300 mm	20/40	24384
15	1100 mm	300 mm	20/40	24385
19	1100 mm	300 mm	20/40	24386
21	1100 mm	300 mm	20/40	24387
24	1100 mm	300 mm	20/40	24388
27	1100 mm	300 mm	20/40	24389




Channel drain

SILICONE FLAT DRAINS

Soft Silicone Flat Drains


Flat drains made of medical grade silicone are used for open and half-open wound drainage by means of gravity and allow a small wound opening thanks to their flat shape.



Product features:

- Soft, medical grade silicone ensures good tissue compatibility
- Fully perforated, with atraumatic openings
- Fully radiopaque
- Flat drain length: 200 mm
- Drain body: 15 Ch., length = 600 mm
- Inner ribbing prevents it from collapsing even at higher pressures or when a vacuum is applied

Flat drain			Drain body			
Width	Height	Length	Ch.	Length	PU/SU	REF
4 mm	3 mm	200 mm	15	600 mm	10/70	21080
7 mm	3 mm	200 mm	15	600 mm	10/70	21081
10 mm	4 mm	200 mm	15	600 mm	10/70	21082
13 mm	5 mm	200 mm	15	600 mm	10/70	21083



Flat drain



ISKUSHEALTH

IMPROVING PATIENT OUTCOMES

Manufacturer

Primed[®]

All rights to technical modifications reserved!

Distributed by:

IRELAND

Iskus Health Ltd

Phone : +353 (0) 1 428 7895 - 96

E-mail: info@iskushealth.com

UNITED KINGDOM

Iskus Health UK Ltd

Phone : +44 (0) 208 869 6509

E-mail: uksales@iskushealth.com

www.iskushealth.com